

Product datasheet for **MG203439**

Nrm (NM_134122) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nrm (NM_134122) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Nrm
Synonyms:	2610307M02Rik; AI429796
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203439 representing NM_134122 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCACCGGCGCTGCTGCTGGTTCCTGCGGCCCTCGCCTCCTTCATCCTGGCCTTTGGCACCGGAGTGG
AATTCGTGCGTTTCACCTCCCTGAGGCCGCTGCTTGGAGGGATCCCGAGTCTGGCGGTCCCGATGCCCC
CCACGGGTGGTTGGCTGCCCTGCAGGACCGGAGCATTCTTGCTCCTTGGCTTGGGATCTATGTCTGCTG
CTACTGTTTGTGGTACAGCACAGCCTTATGGCGACTGAGGCAGTGAAGGCGTGGACATCCCGGTACTTTG
GCGTCCTTCAGAGGTCGCTCTATGTGGCGTGCACAGCCCTGGCCTTGCAATTGGTGATGCGGTACTGGGA
GACTACACCCAGAGGCCAGTGTATGGGAGGCTCGGGCTGAACCATGGGCCACCTGGGTACCTCTCCTC
TGCTTTGTGCTCCACGTTGTTTCCTGGTCCTCATCTTTAGCATCCTTCTGGTCTTTGACTACGTGAAC
TGATGGGCCTCAAGCAGTTTATTACCATGTGCTGGGACTGGGCGAGCCCTGTCCCTGAAGTCTCCTCG
GGCTCTGAGACTCTTTCCACCTTCGACACCCAGTGTGTGGAGCTGCTGACTGTGCTGTGGTGGTC
CCCACCTTGGGCACCGATCGGCTCCTCCTGGCTCTTCTTTACCCTCTACCTGGGCTTGGCTCATGGGC
TTGACCAGCAAGATCTCCGCTACCTTCGGTCTCAGCTACAAAGAAACTCCACCTTCTCTCCAGACCTCA
GGATGGGAAGCAGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >MG203439 representing NM_134122
 Red=Cloning site Green=Tags(s)

MAPALLLVPAALASFILAFGTGVEFVRFTSLRPLLGGIPESGGPDARHGWLAALQDRSILASLAWDLCLL
 LLFVVQHSLMATEAVKAWTSRYFGVLQRSLYVACTALALQLVMRYWETTPRGPVLWEARAEPWATWVPLL
 CFVLHVVSLLIFSILLVFDYAELMGLKQVYYHVLGLGEPLSLKSPRALRLFSLHRHPVCVELLTVLWVV
 PTLGTDRLLLALLFTLYLGLAHGLDQDLRYLRSQQLQRKLHLLSRPQDGEAE

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_134122

ORF Size: 786 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_134122.2](#)

RefSeq Size: 1505 bp


RefSeq ORF: 789 bp

Locus ID: 106582

UniProt ID: [Q8VC65](#)

Cytogenetics: 17 B1

Product images:



Circular map for MG203439