

## Product datasheet for **MG203495**

### Rps3a (BC081451) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Rps3a (BC081451) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rps3a
Synonyms:	Rps3a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203495 representing BC081451 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCGGTCGGCAAGAACAAGCGCCTGACGAAAGGTGGCAAGAAGGGAGCTAAGAAGAAAGTGGTCGATC  
CATTTTCTAAGAAAGACTGGTACGATGTGAAAGCTCCAGCCATGTTCAATATTAGGAACATCGGGAAGAC  
ACTAGTCACGAGGACTCAAGGAACCAAAATTCCTCTGATGGCCTCAAGGGTCGTGTGTTGAAGTGAGC  
CTTGCTGATCTACAGAATGATGAAGTTGCGTTTAGAAAATCAAGCTAATCACTGAGGACGTTCCAGGCA  
AAAAGTGTCTGACTAACTCCATGGTATGGATCTTACCCGGGACAAAATGTGCTCCATGGTCAAGAAGTG  
GCAGACCATGATCGAAGCTCACGTTGATGTCAAGACGACCGATGGGTATTGCTCCGACTTTCTGTGTT  
GGCTTCACTAAGAAACGCAACAACCCAGATAAGAAAGACATCCTATGCGCAGCACCAGCAAGTCCGCCAGA  
TCCGGAAGAAGATGATGGAAATCATGACCCGAGAGGTGCAGACTAATGACTTGAAGGAAGTAGTTAATAA  
ACTGATCCAGACAGCATTGGGAAAGACATAGAAAAGGCTTGCCAGTCCATTTATCCTCTCCATGACGTC  
TTTGTTAGAAAAGTAAAAATGCTGAAGAAACCAAGTTTGAATTGGGCAAGCTCATGGAGCTGCATGGGG  
AAGCGGCAGCTCTGGCAAAGCAGCTGGAGACGAGACGGGTGCTAAGGTTGAGCGAGCTGACGGATATGA  
ACCGCCAGTCCAAGAATCAGTG

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG203495 representing BC081451  
Red=Cloning site Green=Tags(s)

MAVGKNKRLTKGGKKGAKKKVVDPFSSKKDWYDVKAPAMFNIRNIGKTLVTRTQGTKIASDGLKGRVFEVS  
 LADLQNDEVAFRKFKLI TEDVQGNCLTNFHGMDLTRDKMCSMVKKWQTMIEAHVDVKTTDGYLLRFLFCV  
 GFTKKRNNQIRKTSYAQHQQVRQIRKKMMEIMTREVQTNLKEVVKLIPDSIGKDIEKACQSIYPLHDV  
 FVRKVMLKKPKFELGKLMELHGEHGGSSGAAGDETGAKVERADGYEPPVQESV

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shutting:



**Plasmid Map:**



**ACCN:** BC081451

**ORF Size:** 794 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC081451</a> , <a href="#">AAH81451</a>
<b>RefSeq Size:</b>	872 bp
<b>RefSeq ORF:</b>	794 bp
<b>Locus ID:</b>	20091
<b>Cytogenetics:</b>	3 F1
<b>Gene Summary:</b>	May play a role during erythropoiesis through regulation of transcription factor DDIT3. [UniProtKB/Swiss-Prot Function]