

Product datasheet for **MG203535**

Rpl7a (BC084678) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rpl7a (BC084678) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rpl7a
Synonyms:	Surf3, Surf-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203535 representing BC084678 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCCAAAGGGAAGAAGGCCAAGGGGAAGAAGGTGGCCCCGGCCCCCGCGTCGTCAAGAAGCAGGAGG
CCAAGAAGGTGGTCAACCCTCTGTTTCGAGAAGAGGCCCAAGAAGCTTCGGCATTGGACAGGACATCCAGCC
CAAAAGAGATTTAACGCGCTTCGTCAAATGGCCCCGCTACATCAGACTGCAGCGGCAAAGAGCTATCCTC
TACAAGCGGCTCAAAGTCCCTCCTGCCATTAACCAGTTACCCAGGCCCTGGACAGGCAGACAGCTACTC
AGCTGCTTAAGCTTGCCACAAGTACAGGCCAGAGACAAAGCAAGAGAAGAAGCAAGGCTACTGGCCCG
TGCTGAGAAGAAAGCTGCTGGCAAAGGGGACGTCCTCAACTAAGAGACCACCTGTCCTCCGAGCAGGAGTC
AATACAGTCAACCCTTGGTGGAGAACAAGAAGGCTCAGCTGGTGGTATTGCCCATGACGTAGACCCCA
TTGAGCTGGTGGTTTTCTACCTGCTCTGTGTCGAAAGATGGGGGTGCCCTACTGCATCATCAAGGGAAA
GGCCAGGCTGGGGCACCTGGTCCACAGGAAGACATGCACCACCGTTGCCCTCACACAGGTTAACTCGGAA
GACAAGGTGCTCTGGCTAAGCTGGTGGAAAGCTATTAGGACCAATTATAATGACAGATATGACGAGATCC
GTCGCCACTGGGGAGGCAACGTCCTGGTCCCTAAGTCTGTGGCTCGAATTGCCAAGCTGGAAAAAGCAA
GGCTAAAGAACTCGCCACTAAATTGGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203535 representing BC084678
 Red=Cloning site Green=Tags(s)

MPKGGKAKGKKVAPAPAVVKKQEAKKVVNPLFEKRPKNFGIGQDIQPKRDLTRFVKWPRYIRLQRQRAIL
 YKRLKVPPAINQFTQALDRQATATQLLKLAKHYRPETKQEKQRLARA EKKAAGKGDVPTKRPPVLRAGV
 NTVTTLVENKKAQLVVI AHDVDP IELVVFLPALCRKMGVPYCI IKGKARLGH LVHRKTCTTVAFTQVNSE
 DKGALAKLVEAIRTNYNDRYDEIRRHWGNNVLGPKSVARIAKLEKAKAKELATKLG

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC084678

ORF Size: 800 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC084678 , AAH84678
RefSeq Size:	897 bp
RefSeq ORF:	800 bp
Locus ID:	27176
Cytogenetics:	2 19.1 cM