

Product datasheet for **MG203592**

Mxi1 (BC064453) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mxi1 (BC064453) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mxi1
Synonyms:	bHLHc11; ENSMUSG00000067085; Gm10197; Mad2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203592 representing BC064453 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCCCCGCCGCGCCGCGCCCCAGCCCCGGCGCAGCCGGAGGAGCCGGCGGGGGCCAAGCCCCGGT
GCCCTTCTCGGACATTTTCAACACCAGCGAGAAGCTCGATGGAGAAGCACATCAACACTTTTCTGCAGAA
CGTGCAGATTCTGCTCGAGGCAGCCAGCTACCTGGAGCAGATCGAGAAAGAAAACAAAAAGTGTGAACAT
GGCTACGCCTCATCGTTCCCTCCATGCCGAGCCCCGGCTACAGCACTCGAAGCCCCACGGAGGTTGA
GCCGGGCACAGAAACACAGCAGTGAAGCAGCAACACCAGCACTGCCAACAGATCTACACACAATGAGTT
GGAAAAGAACCGACGAGCTCACCTGCGCCTGTGTTTAGAACGCTTGAAAGTTCTGATCCCGCTGGGCCCA
GACTGCACCAGGCACACAACACTCGGTTTGTCTCAACAAAGCCAAAGCACACATCAAGAACTTGAAGAAG
CGGAGAGGAAGAGCCAGCACCAGCTAGAGAAGTTGGAACGAGAACAGAGGTTTTTAAAGCGCGACTGGA
ACAGCTGCAGGGGCTCAGGAGATGGAGCGGATACGAATGGACAGCATTGGATCAACCATCTCTTCAGAT
CGCTCGGATTCAGAGCGAGAGGAGATTGAAGTGGATGTGGAAAGCACAGAGTTCTCCCATGGAGAAGCAG
ACAGTGTGAGTACCACCAGCATCAGTGACCTTGACGACCACAGCAGCCTGCAGAGTGTGGGAGTGACGA
GGGTTATTCCAGTGCCAGTGTCAAACCTCTCTTCGCGTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203592 representing BC064453
Red=Cloning site Green=Tags(s)

MPPAAAAPQPPAQPEEPAGAKPRCPFSDIFNTSENSMEKHINTFLQNVQILLEAASYLEQIEKENKKCEH
 GYASSFSPMSPRLQHSKPPRRLSRAQKHSSGSSNTSTANRSTHNELEKNRRAHLRLCLERLKVLIPLGP
 DCTRHTTLGLLNKAKAHIKKLEEAERKSQHQLENLEREQRFLLKRRLEQLQGPQEMERIMDSIGSTISSD
 RSDSEREEIEVDVESTEFSHGEADSVTTSISDLDDHSSLQSVGSDEGYSSASVKLSFAS

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC064453

ORF Size: 812 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC064453 , AAH64453
RefSeq Size:	2643 bp
RefSeq ORF:	812 bp
Locus ID:	17859
Cytogenetics:	19 47.53 cM
Gene Summary:	This gene encodes a protein containing a helix-loop-helix domain characteristic of transcription factors, which allows heterodimerization and sequence-specific DNA binding. The encoded protein is related to a family of Myc/Max/Mad proteins that are involved in the regulation of several cellular processes. The protein encoded by this gene is a transcriptional repressor thought to negatively regulate Myc function. Three alternatively spliced transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]