

## Product datasheet for **MG203625**

### Cela2a (NM\_007919) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Cela2a (NM\_007919) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Cela2a  
**Synonyms:** El; Ela; Ela-2; Ela2; Ela2a  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG203625 representing NM\_007919  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATCAGGACACTGCTGCTATCTGCCTTGGTGGCTGGAGCCCTCAGCTGTGGGTACCCACTTATGAGG  
TGGAGGATGATGTGAGCAGGGTAGTTGGGGTCAAGAGGCCACACCAACACCTGGCCCTGGCAGGTCTC  
CCTGCAGGTCCTTTCCTCCGGAAGGTGGCGCCACAACCTGCGGAGGCTCCCTGGTGCCAACTGGGTT  
CTGACAGCTGCCATTGCCTCAGCAACTATCAGACCTACCGAGTGTCTGGGCGCACACAGCCTCTCCA  
ACCCCGAGCTGGCTCTGCTGCTGTTCAAGTCTCTAAGCTTGTGGTCCACCAGAGGTGGAACCTCCAAAA  
CGTCGGCAATGGCTATGACATTGCCTTAATCAAACCTGGCCAGCCAGTGACCCTGAGCAAGAATCCAG  
ACAGCTTGCTCCACCCGCTGGCACCATTCTCCGAGAACTATGTCTGCTATGTCACAGGCTGGGGCC  
TGCTGCAGACCAATGGGAACAGTCTGACACCCTGAGGCAGGGCCGCTGCTGGTGTGGACTATGCCAC  
CTGCTCCAGCGCTAGCTGGTGGGAAAGCTCTGTGAAGTCCAGCATGGTGTGCGCTGGTGGCAGCGCGT  
ACCTCCAGCTGCAATGGGACTCTGGCGGACCACTGAATTGCCGGCATCTAATGGCCAGTGGCAGGTGC  
ATGGCATCGTGAGCTTCGGCTCCTCTCTGGGCTGCAACTACCCCGCAAGCCATCCGCTTCACCAAGGT  
CTCCAACATGACTGGATCAACTCGGTGATGGCAAGGAAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG203625 representing NM\_007919  
 Red=Cloning site Green=Tags(s)

MIRTL L L S A L V A G A L S C G Y P T Y E V E D D V S R V V G G Q E A T P N T W P W Q V S L Q V L S S G R W R H N C G G S L V A N N W  
 L T A A H C L S N Y Q T Y R V L L G A H S L S N P G A G S A A V Q V S K L V V H Q R W N S Q N V G N G Y D I A L I K L A S P V T L S K N I Q  
 T A C L P P A G T I L P R N Y V C Y V T G W G L L Q T N G N S P D T L R Q G R L L V D Y A T C S S A S W W G S S V K S S M V C A G G D G V  
 T S S C N G D S G G P L N C R A S N G Q W Q V H G I V S F G S S L G C N Y P R K P S V F T R V S N Y I D W I N S V M A R N

TRTRPLE - GFP Tag - V

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja2368\\_g09.zip](https://cdn.origene.com/chromatograms/ja2368_g09.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_007919


**ORF Size:** 813 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_007919.3</u>
<b>RefSeq Size:</b>	912 bp
<b>RefSeq ORF:</b>	816 bp
<b>Locus ID:</b>	13706
<b>UniProt ID:</b>	<u>P05208</u>
<b>Cytogenetics:</b>	4 74.75 cM
<b>Gene Summary:</b>	This gene encodes a serine protease enzyme that hydrolyzes elastin. This gene is highly expressed in the pancreatic acinar cells where the encoded preproprotein undergoes processing including signal peptide cleavage to generate an inactive zymogen. The removal of N-terminal activation peptide from the zymogen by trypsin generates active elastase enzyme. This gene is also expressed in the mouse epidermis where it participates in pro-filaggrin processing. [provided by RefSeq, Jul 2016]

**Product images:**


Circular map for MG203625