

Product datasheet for **MG203695**

Odf3l1 (BC049697) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Odf3l1 (BC049697) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Odf3l1
Synonyms:	MGC58606
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203695 representing BC049697 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCAGTCCAAGGGTCCAGGAACATGTGTTCTACGCACAGCACCCAGAGAAGGAGGAAGAGGTGC
CCTCATGGCATGAAATCAAGCAGACCCCTGTATCATGGCCACAATCAAAGGTCCAGGGCCTGCCAAGTA
CCTCCGGTCATCCTGTACTGGCTATATAGCCCATGACATCTCCATGTTCCAGGAGCCAGCTTACAGTCTG
CATACACGGCACACTAAGAAACGGATCATAGACATAAACAGCCCCGACCTTGCTATTTCTTGAATCCCA
AAGTAACTCGTTTGGAAATATCCACCTGCCCCAGGTCCCATGGAGGAGCGCATCTCCAATCCACGCAT
AAATTGCATGCCAGCCTCCTGCAAGTACAACCTTGAGAAGACTCGACCCTCTGGGGAACGTGAGCCTCCC
CAGTACACTTTTGGCTACCGGTGTCATATAGAGTGATGGATCCCAACCTGCTCCTAACCAAGTACCAGC
TGCCAGTCACTCTGGGACCAACATCCCTGTCTTCCGAGCTGCCCTTCTATAGTCTGGCTTCCACAAA
CAAGAATTGGTTCACAAGGAGAATATAGCAGGAGGTCCCGACCTGCTATGCATACCCGACCTGAGCCA
TCCGTCTACCAGAACCGAGCCCTTTGTTAGCATGGCCAAGCGCTTTGGCTGTCCTTTGGACCACACGC
ATCGGCCTGGCCCTGGCTCCCACGACATCCAGCCGGTCACTGTGCACAAACCCGCATACCTGCCTTCAC
CATGGGCATCAAACATTCACCCCACTTGTGCCCACTGATTGTGGACATTTGTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >MG203695 representing BC049697
 Red=Cloning site Green=Tags(s)

MKQSKGSRNYVFYAQHPEKEEEVPSWHEIKQTPVIMATIKGPGPAKYLRSSCTGYIAHDISMFQEPAYSL
 HTRHTKKRIIDINSPGPCYFLNPKVTRFGISTCPQVPMEEERISNPRINCMASCKYNLEKTRPSGERQPP
 QYTFGYRCPYRMDPNPAPNQYQLPVTLTGNIPVFRAAPSYSLASTNKNWFHKENIAGGPGPAMHTRPEP
 SVYQNRSPFLFMAKRFGCPLDHTHRPGPGSHDIQPVTVHKPRIPAFTMGIKHSPHLCPLIVDIDC

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC049697

ORF Size: 827 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: BC049697, AAH49697


RefSeq Size: 1204 bp

RefSeq ORF: 827 bp

Locus ID: 382075

Cytogenetics: 9 B

Product images:



Circular map for MG203695