

Product datasheet for **MG203701**

Tpsb2 (NM_010781) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tpsb2 (NM_010781) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Tpsb2
Synonyms:	AV011504; Mcp-6; Mcp6; Mcpt6; MMCP-6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203701 representing NM_010781 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTGAAGCGGGCTGCTGCTGTGGCACTGTCCCTCCTGGCTAGTCTGGTGTACTCAGCCCTC
GCCAGCCAATCAGCGAGTGGGCATCGTGGGAGGACATGAGGCTTCTGAGAGTAAGTGGCCCTGGCAGGT
GAGCCTGAGATTTAAATTAACACTGGATACATTTCTGCGGAGGCTCTCTCATCCACCCACAGTGGGTG
CTCACTGCGGCACACTGTGTGGACCGCACATCAAAGCCACAGCTCTCCGGGTGCAGCTTCGTGAGC
AGTATCTATACTATGGGACCAGCTCCTCTCTTGAACCGGATCGTGGTGCACCCCACTATTACACGGC
CGAGGGTGGGGCAGACGTTGCCCTGCTGGAGCTTGAGGTCCTGTGAATGTCTCCACCCATATCCACCC
ATATCCCTGCCCCCTGCCTCGGAGACCTTCCCCCTGGGACATCGTGTGGGTGACAGGCTGGGGCGACA
TTGATAATGACGAGCCTCTCCACCTCCTATCCTCTGAAGCAAGTGAAGTTCCCATTTGTGAAAAACAG
CCTGTGTGACCGAAGTACCACACTGGCCTCTACACGGGAGATGATTTCCCATTTGTCATGATGGCATG
CTGTGTGCTGAAATACCAGGAGAGACTCCTGCCAGGGCGATTGAGGGGGCCACTGGTCTGCAAAGTGA
AGGGTACCTGGCTGCAGGCAGGAGTGGTCAGCTGGGGTGGGGTGCAGCAGCCCAACAAGCCTGGCAT
CTACACCCGGGTGACATACTACTTAGACTGGATCCACCGCTATGTCCCTGAGCATTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203701 representing NM_010781
 Red=Cloning site Green=Tags(s)


MLKRLLLLLWALLSLLASLVYSAPRPNQRVGVGGHEASESKWPWQVSLRFKLNWYIHFCCGSLIHPQWV
 LTAHCVGPHIKSPQLFRVQLREQYLYGDQLLSLNRIVVHPHYTAEGGADVALLELVPVNVSTHIHP
 ISLPPASETFPPGTSCWVTGWGDIDNDEPLPPPYPLKQVKVPIVENS�CDRKYHTGLYTGDDFPIVHDGM
 LCAGNTRRDSCQGDSSGGPLVCKVKGTWLQAGVVSWGEGCAQPNKPGIYTRVITYLLDWIHRVYVPEHS

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_010781

ORF Size: 828 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010781.2
RefSeq Size:	1088 bp
RefSeq ORF:	831 bp
Locus ID:	17229
Cytogenetics:	17 12.53 cM