

Product datasheet for **MG203704**

Ctdspl (BC094289) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ctdspl (BC094289) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ctdspl
Synonyms:	SCP3, HYA22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203704 representing BC094289 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGGTCCGGCCATCATCACCCAGGTGACCAACCCCAAGGAGGATGAGGCCCGTCCGCCGTCGCGG
GCGAAAAAGCTTCCAGCGCAACATCAGCCTAAAGAAGCAAAGGGGTCGGAGCATCCTTAGCTCTTTCTT
CTGCTGCTCCGAGACTACAATGTGAAGCCCCTCCGGCCAACAGCCCCAGTGTGCTCCGCCACTGGTG
GAGGAGAACGGTGGGCTCAGAAGGGTGACCAGAGGCAGGTCAATCCCGTACCAAGCCCACCAGCTAAAT
ACCTCCTCCAGAGGTGACGGTCTCGACTATGGAAAGAAATGTGTGGTCATTGATCTAGACGAAACCT
GGTGACAGTTCGTTAAGCCCATTAGTAATGCTGACTTCATTGTTCTGTTGAAATCGACGGAACCATA
CACCAGGTGTATGTGCTGAAGCGACCTCAGTGGACGAGTTCCTCCAGAGGATGGGGCAGCTCTTTGAGT
GTGTGCTCTTTACTGCCAGCTTGCCCAAGTATGCAGACCTGTTGCCGACCTCCTGGACCCTGGGGTGT
GTTCCGGGCCGACTCTCAGAGAATCCTGCGTGTCCACCGTGGGAACTATGTGAAGGACTTGAGTCGC
CTGGGACGGGAGCTGAGCAAAGTATCATCGTGGACAATTCTCCTGCCTCCTACATCTTCCATCCTGAGA
ATGCAGTGCCCGTGCAGTCTGGTTTGACGACATGACAGACCGAGCTGTTGGACCTCATCCCGTTCTT
CGAGGGCTGAGCAGGGAAGACGATGTCTACAGCATGCTGCACAGACTCTGTAGCAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203704 representing BC094289
 Red=Cloning site Green=Tags(s)

MDGPAIITQVTNPKEDARSVPVAGEKASQRNISLKKQRGRSILSSFFCCFRDYNVEAPPANSPSVLPPLV
 EENGLQKGDQRQVIPVPSPPAKYLLPEVTVLVDYGKCCVVIDLDELTVHSSFKPISNADFIVPVEIDGTI
 HQVYVLRPHVDEFLQRMGQLFECVLFATSLAKYADPVADLLDRWGVFRARLFRESCVFRHNGYVKDLRSR
 LGRELSKVIIVDNSPASYIFHPENAVPVQSWFDDMTDELLDLIPFFEGLSREDDVYSMLHRLCSR

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: BC094289

ORF Size: 830 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC094289 , AAH94289
RefSeq Size:	1051 bp
RefSeq ORF:	830 bp
Locus ID:	69274
Cytogenetics:	9 F3
Gene Summary:	Preferentially catalyzes the dephosphorylation of 'Ser-5' within the tandem 7 residue repeats in the C-terminal domain (CTD) of the largest RNA polymerase II subunit POLR2A. Negatively regulates RNA polymerase II transcription, possibly by controlling the transition from initiation/capping to processive transcript elongation. Recruited by REST to neuronal genes that contain RE-1 elements, leading to neuronal gene silencing in non-neuronal cells (By similarity).[UniProtKB/Swiss-Prot Function]