

## Product datasheet for **MG203850**

### Six1 (NM\_009189) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Six1 (NM\_009189) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Six1  
**Synonyms:** BB138287  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG203850 representing NM\_009189  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGTCGATGCTGCCGTCGTTTGGTTTTACGCAAGAGCAAGTGGCGTGCCTGTGCGAAGTTCTGCAGCAAG  
GAGGGAACCTGGAACGCCTGGGCAGGTTCTTGTGGTCGTTGCCCGCCTGCGATCACCTGCACAAGAACGA  
GAGCGTGCTCAAGGCCAAGGCGGTGGTCGCCTCCACCGCGCAACTCCGCGAGCTCTACAAGATACTG  
GAGAGCCACCAGTTCTCGCTCACAATCACCCAAACTGCAGCAGCTGTGGCTGAAAGCGCACTACGTGG  
AGGCCGAGAACTTCGCGCCGACCCCTGGGTGCCGTGGGCAAAATATCGGGTGCGCCGAAAATCCCGTT  
GCCCGGACCATCTGGGACGGCGAGGAGACCAGCTACTGCTTTAAGGAGAAGTCTCGGGGCGTGCTGCGG  
GAGTGGTACGCGCACAACCCCTACCCCTCACCGAGGGAGAAACGGGAGCTGGCCGAGGCCACCGGCCTCA  
CCACCACCAGGTCAGCAACTGGTTTAAGAACCGAGGCAAGAGACCGGGCCGCGAGGCCAAGGAAAAG  
GGAGAACACCGAAAACAATAACTCCTCCTCCAACAAGCAGAATCAACTCTCCTCTGGAAGGGGCAAG  
CCGCTCATGTCCAGCTCAGAAGAGGAGTTCTCACCCCCAAAGTCCAGACCAGAACTCGGTCTTCTGC  
TCCAGAGCAATATGGGCCACGCCAGGAGCTCAAATTTCTTCCAGGCCTCACAGCCTCCAGCCAG  
CCACGGTCTGAAGCCCATCAGCACCAGCTCCAGGACTCTGCTGGGCCACTCACCTCCAGTTTGGTG  
GACTTGGGTTCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203850 representing NM\_009189  
 Red=Cloning site Green=Tags(s)

MSMLPSFGFTQEQVACVCEVLQQGGNLERLGRFLWSPACDHLHKNESVLKAKAVVAFHRGNFRELKIL  
 ESHQFSPHNHPKLQQLWKAHYVEAEKLRGRPLGAVGKYRVRKFPPLRTIWDGEETSYCFKEKSRGVL  
 R E WYAHNPYPSPREKRELAETGLTTTQVSNWFKNRRQRDRAAEAKERENTENNSSSNKQNLSPLEGGK  
 PLMSSSEEEFSPQSPDQNSVLLQLQSNMGHARSSNYSLPGLTASQPSHGLQAHQHQLQDSSLGPLTSSLV  
 DLGS

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM\_009189

ORF Size: 852 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_009189.3](#)

**RefSeq Size:** 2601 bp

**RefSeq ORF:** 855 bp

**Locus ID:** 20471

**UniProt ID:** [Q62231](#)

**Cytogenetics:** 12 30.34 cM

**Gene Summary:**

Transcription factor that is involved in the regulation of cell proliferation, apoptosis and embryonic development. Plays an important role in the development of several organs, including kidney, muscle and inner ear. Depending on context, functions as transcriptional repressor or activator. Lacks an activation domain, and requires interaction with EYA family members for transcription activation. Mediates nuclear translocation of EYA1 and EYA2. Binds the 5'-TCA[AG][AG]TTNC-3' motif present in the MEF3 element in the MYOG promoter. Regulates the expression of numerous genes, including MYC, CCNA1, CCND1 and EZR. Acts as activator of the IGFBP5 promoter, probably coactivated by EYA2. Repression of precursor cell proliferation in myoblasts is switched to activation through recruitment of EYA3 to the SIX1-DACH1 complex. During myogenesis, seems to act together with EYA2 and DACH2. [UniProtKB/Swiss-Prot Function]