

Product datasheet for **MG203869**

Plet1 (NM_029639) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Plet1 (NM_029639) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Plet1
Synonyms:	0610037B23Rik; 1600029D21Rik; AA407923; AU043048
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203869 representing NM_029639 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGCTGTCGCTCCGCTCCTTGCTTCCACACCTGGGACTGTTCTGTGCCTGGCTCTGCACTTATCCCCCT
CCCTCTCTGCCAGTGATAATGGTCCTGCGTGGTCCTTGATAACATCTACACCTCCGACATCTTGGAAAT
CAGCACTATGGCTAACGTCTCTGGTGGGATGTAACCTATACAGTGACGGTCCCCGTGAACGATTTCAGTC
AGTGCCGTGATCCTGAAAGCAGTGAAGGAGGACGACAGCCAGTGGGCACCTGGAGTGGAACATATGAGA
AGTGCAACGACAGCAGTGTCTACTATAACTTGACATCCCAAAGCCAGTCGGTCTTCCAGACAACTGGAC
AGTTCTACTTCCGAGGATGTGACTAAAGTCAACCTGCAGGTCCCTCATCGTCGTAATCGCACAGCCTCA
AAGTCATCCGTGAAAATGGAACAAGTACAACCCCTCAGCCTCAACCCCTATTCTGAGAGTTCTGAGACCA
GCCAGACCATAAACACGACTCCAAGTGAACACAGCCAAGACTACAGCCAAGGACACAGCCAATACCCAC
AGCCGTGACCACAGCCAATACCACAGCCAATACCACAGCCGTGACCACAGCCAAGACCACAGCCAAAAGC
CTGGCCATCCGCACTCTCGGCAGCCCCCTGGCAGGTGCCCTCCATATCCTGCTTGTCTTCTCATTAGTA
AACTCCTCTTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >MG203869 representing NM_029639
 Red=Cloning site Green=Tags(s)

MLSLRSLPHLGLFLCLALHLSPSLASDNGSCVLDNIYTSIDILEISTMANVSGGDVITYTVTPVNDVS
 SAVILKAVKEDDSPVGTWSGTYEKCNDSVYYNLTSSQSVFQTNWTVPTSEDVTKVNLQVLIVNRTAS
 KSSVKMEQVQPSASTPIPESETSQITINTPTVNTAKTTAKDTANTTAVTTANTTANTTAVTTAKTTAKS
 LAIRTLGSPLAGALHILLVFLISKLLF

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_029639

ORF Size: 711 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_029639.2](#), [NP_083915.2](#)

RefSeq Size: 1846 bp

RefSeq ORF: 714 bp


Locus ID: 76509

UniProt ID: [Q8VEN2](#)

Cytogenetics: 9 A5.3

Gene Summary: Modulates leading keratinocyte migration and cellular adhesion to matrix proteins during a wound-healing response and promotes wound repair. May play a role during trichilemmal differentiation of the hair follicle.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG203869