

## Product datasheet for **MG203904**

### Stx3 (BC024844) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Stx3 (BC024844) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Stx3
Synonyms:	Syn-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203904 representing BC024844 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAAGGACCGGCTGGAGCAGCTGAAGGCCAAGCAGCTGACGCAGGATGATGACACGGACGAGGTTGAGA  
TTGCCATTGACAACACAGCCTTCATGGACGAGTTCTTCTCTGAGATCGAGGAAACACGGCTCAACATTGA  
CAAGATCTCGGAGCATGTAGAGGAAGCTAAGAACTCTACAGTATCATTCTCTGCACCGATCCCAGAG  
CCAAAAACCAAAGATGACCTTGAACAGCTCACAACCTGAGATCAAGAAAAGGGCCAACAACGTCCGGAACA  
AGCTGAAGAGCATGGAGAAGCACATTGAGGAAGATGAGGTTCCGGTCATCTGCAGACCTTCGGATACGAAA  
GTCCCAGCACTCCGTCCTCTCCCGAAGTTTGTGGAGGTGATGACAAAGTACAATGAGGCTCAGGTGGAC  
TTCCGTGAACGCAGCAAAGGGCGCATCCAACGGCAGCTTGAAATTAAGTGGCAAGAAGACAACAGATGAGG  
AGCTGGAAGAGATGTTGGAGAGTGGCAACCCAGCCATCTTCACTTCTGGGATCATTGACTCCCAGATTTT  
CAAGCAAGCCCTCAGCGAGATTGAGGGTCGGCACAAGGACATCGTAAGGCTGGAGAGCAGCATCAAGGAG  
CTCCATGACATGTTTATGGACATCGCCATGCTGGTGGAAAATCAGGGAGCCATGATTGACCGCATTGAGA  
ACAACATGGACCAGTCAGTGGGCTTTGTGGAACGGGCTGTGGCAGATACCAAAAAGGCGGTCAAGTATCA  
GAGTGAAGCCCGGAGGAAGAAGATCATGATCATGATCTGCTGTATTATCCTTGCCATCATCTTGGCTTCT  
ACCATTGGGGGCATATTTGCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG203904 representing BC024844  
Red=Cloning site Green=Tags(s)


MKDRLEQLKAKQLTQDDTDEVEIAIDNTAFMDEFFSEIEETRLNIDKISEHVVEAKKLYSIILSAPIPE  
 PKTKDDLEQLTTEIKKRANVRNKLKSMEKHIEEDEVSRSSADLRIRKSQHSVLSRKFVEVMTKYNEAQVD  
 FRERSKGRIQRQLEITGKKTDEELEEMLESGNPAIFTSGIIDSQISKQALSEIEGRHKDIVRLESSIKE  
 LHDMFMDIAMLVENQGMIDRIENNMDQSVGFVERAVADTKKAVKYQSEARRKKIMIMICCIILAILAS  
 TIGGIFA

TRTRPLE - GFP Tag - V


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** BC024844

**ORF Size:** 863 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC024844</a> , <a href="#">AAH24844</a>
<b>RefSeq Size:</b>	2495 bp
<b>RefSeq ORF:</b>	863 bp
<b>Locus ID:</b>	20908
<b>Cytogenetics:</b>	19 A
<b>Gene Summary:</b>	Potentially involved in docking of synaptic vesicles at presynaptic active zones. [UniProtKB/Swiss-Prot Function]