

Product datasheet for **MG203921**

Dhrs7 (NM_025522) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dhrs7 (NM_025522) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Dhrs7
Synonyms:	2310016E22Rik; 5730564L20Rik; AW061210; Retdsr4; Retsdr4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203921 representing NM_025522 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTGGTGTGGGTTACTGGAGCGTCAAGTGGCATTGGCGAGGAGCTAGCTTCCAGTTATCTAACTCG
GGGTCTCTGGTGTCTCAGCCGAAGGGCGCAGGAGCTGGAGCGGGTGAAGAGAAGATGCCTGGAAAA
TGGCAATTTAAAAGAAAAGATACTGGTTCTGCCCTTGACTTGACCGATAACAAGTCCCATGAGGCA
GCCACCAAGCCGTTCTGCAGGAGTTTGGTAAAATTGACATTTGGTCAACAACGGTGAAGAAGCCAGC
GTTCCCTTAGTTCTAGAAACCACTTGGATGTCTTCAAGGAGCTGATAAATCTGAACATACATAGGGACAGT
ATCCTTGACCAAGTGTGTCTGCCTCACATGATAGAGAGGAAACAAGGGAAGATTGTGACTGTGAACAGC
ATCGCGGGGATTGCCTCTGTGTCCCTTCCAGTGGCTACTGTGCCAGCAAGCATGCTCTTCGGGGTTTTT
TCAATGCCTCCATAGTGAGCTCGGCCAATATCCAGGCATAACATTTTGCAACGTTTACCCGGGGCCTGT
GCAGTCCGACATTGTGAAGAATGCCTTACCGAGGAAGTAACTAAGTCTATGAGAAACAACATAGACCAG
AGTTACAAGATGCCAACGAGCCGCTGTGTGCGGCTGATGCTCATCTCCATGGCCAATGATCTGAAGGAAG
TCTGGATCTTGACCATCCTGTGCTGTTGGGGCATACATATGGCAGTATATGCCAACCTGGGCCGATG
GCTAAACTGCAAGTTGGGAAGGAAAGAATCCAAAACCTTAAGAATAATTTGGATCCAGACTTGCCTTAT
AAATTCCTGAAGCAAAGAAGGAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203921 representing NM_025522
 Red=Cloning site Green=Tags(s)

MVVVWTGASSGIGEELAFQLSKLGVSLVLSARRAQELERVKRRCLEENGLKEKDILVPLDLTDTSSHEA
 ATKAVLQEFGKIDILVNNGGRSQRSLVLETNLDVFKELINLNYIGTVSLTKCVLPHMIERKQGKIVTVNS
 IAGIASVLSSSGYCASKHALRGFFNALHSELGQYPGITFCNVYPGPVQSDIVKNAFTEEVTKSMRNNIDQ
 SYKMPTSRCVRLMLISMANDLKEVWISDHPVLLGAYIWQYMPYWAAWLNCKLGKERIQNFKNLDPDLPLP
 KFLKAKKD

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_025522

ORF Size: 864 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_025522.4
RefSeq Size:	1256 bp
RefSeq ORF:	1017 bp
Locus ID:	66375
UniProt ID:	Q9CXR1
Cytogenetics:	12 C3