

Product datasheet for **MG203952**

C1qtnf7 (NM_175425) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C1qtnf7 (NM_175425) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	C1qtnf7
Synonyms:	5530401N20Rik; 8430425G24Rik; Ctrp7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG203952 representing NM_175425 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATTGCTCTACGTGACGAGTCTGCCATCTGTGCAAGTGGACAACCTCGGGCCAATCAGGCTA
AGGGAGAGAGCTACTCTCCAAGGTACATCTGCAGCATCCCTGGATTACCTGGCCCCCAGGTCCTCTGG
AGCAATGGCTCCCCTGGACCCCATGGTCGCATTGGCCTTCTGGAAGGGATGGTAGAGATGGCAGAAAA
GGAGAGAAGGGGAAAAGGGCACTGCAGGTCTAAAAGGTAAAACCTGGACCCTGGCCTTGCTGGTGAGA
AAGGAGACCAAGGAGAACTGGAAGAAAGGACCCATAGGACCAGAGGGTGAAGGAGAAAGTCCGGTCC
AGCTGGGCTCCTGGGCAAAGGGAGACAGAGGAGATCAAGGGGACCCAGGGCTGCCTGGAGTGTGCAGG
TGTGGAAGCATTGTGCTCAAATCTGCCTTTTCAGTTGGCATCACAAACAGCTACCCAGAAGAAAGACTAC
CCATCATATTTAACAAGTCTCTTCAATGAGGGGGAGCATTACAACCTGCAACGGGGAAGTTCATTTG
CGCTTTCCAGGGATCTATTACTTTTCTTATGACATCACGTTGGCCAATAAGCACCTAGCAATCGGGCTG
GTGCACAATGGGCAGTACCGGATAAGGACCTTTGATGCCAACACAGGGAACCATGATGTGGCGTCGGGT
CCACAGTCATCTACCTGCAGCCAGAAGATGAGGTCTGGCTGGAGATCTTCTTCAATGACCAGAACGGCCT
CTTCTCAGATCCAGGCTGGCAGACAGCTTGTCTCTGGGTTTCTCTCTATGTCGATACAGATTACCTG
GATTCTATATCAGAGGATGATGAGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG203952 representing NM_175425
Red=Cloning site Green=Tags(s)

MIVLLYVTSLAICASGQPRANQAKGESYSPRYICSIPLPGPPGPPGANGSPGPHGRIGLPGRDGRDGRK
 GEKGEKGTAGLKGKTGPLGLAGEKGDQGETGKKGPIGPEGEKGEVGPAGPPGPKGDRGDQGDPLPGVCR
 CGSIVLKSAFVSGITTSYPEERLPPIIFNKVLFNEGEHYNPATGKFCAFPGIYYFSYDITLANKHLAIGL
 VHNGQYRIRTFDANTGNHVDVASGSTVIYLPQPEDEVWLEIFFNDQNLGFLSDPGWADSLFSGFLLYVDTDYL
 DSISEDDEL

TRTRPLE - GFP Tag - V

Restriction Sites:


SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_175425

ORF Size: 1292 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_175425.4 , NP_780634.2
RefSeq Size:	3843 bp
RefSeq ORF:	870 bp
Locus ID:	109323
UniProt ID:	Q8BVD7
Cytogenetics:	5 B3