

Product datasheet for **MG204003**

Mrpl19 (NM_026490) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrpl19 (NM_026490) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mrpl19
Synonyms:	9030416F12Rik; 9130412E02Rik; D6Ertd157e; MRP-L15; RLX1; Rpml15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204003 representing NM_026490 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCCTCCATGGCAGAAAGTTGTAGGGCGTCCCTGTACCTGGCCCGGAGCGTCCGAATGGCCAGGC
CACGGCTTGTGCCTTCGCATCTGATGCCTGCCGTGTCTGCACGGGCCAAGCCGATTTCCAGAGCACCGG
GCCTTCAGAGCCCGGTGGGTTCAAGCCACCGCCGAAACCTGTATCGTGGACCGGCCGAGTCCCGGAA
GACGAGAGAAGGTTCTGAGTCTGAATTCATCCCTCCAGGGGGAACAAACCCTCTGAAATTTAAGA
TAGAAAGAAAAGATATGCTAGACAGGAGAAAGGTGCTTCCGATTCAGAGTTCTATGTTGGAAGCATACT
TCGAGTCAACCACAGCCGACCCCTATGCCAGTGGAAAAACCGCCAGTTCCTGGGGATTTGCATCAAGCGA
TCAGGGAATGGGCTCGGGGCTACTTTCCACCTCAGGAACACGATTGAAGGACAAGGTGTGGAGATTTGCT
TTGAACTGTATAATCCTCGAATCCAGGAGATCCAGGTGGTGAAGTTGGAGAAGCAGCTGGATGACAACT
GCTGTACTTACGGGACGCCCTTCCCGAGTACAGCACCTTTGACGTGAACATGAAGCCGGTCCCACAGGAG
GCCTGCCAGGAAGTGCCTGTTAAACAAGCTGAAAGTTAAATGAAGCCAAAGCCGTGGTCCAAACGCTGGG
AACGTCCAAACTCAACATTAAGGGAATCAGATTTGACCTTGCTTAACTGAAGACAAAATGAAGGAAGC
CCAGAAGTGGAACAAGCCATGGATTGAGTTTGATGATGAGAGAATATGACACTTCAAAAATTTGAAGCT
GCATTGTGGGAAGAAATAGAAGCATCGAAGAAGTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG204003 representing NM_026490
Red=Cloning site Green=Tags(s)

MAASMAESCRASLYLARSVRMARPRLAASFADACRVCTGPSRFQSTGPSEPGGFKPPPKPVIIVDRRRVPE
 DERRFLSPEFIPPRGRNPLKFKIERKDMLDRRKVLPIPEFYVGSILRVTTADPYASGKTSQFLGICIKR
 SGNLGLATFTLRNTIEGQVEICFELYNPRIQEIQVVKLEKQLDDNLLYLRLDALPEYSTFDVNMKPVPQE
 ACQEVVNNKLVKMKPKPWSKRWERPNFNKIGIRFDLALTEEQMKEAQKWNKPWIEFDMMREYDTSKIEA
 ALWEEIEASKKS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_026490

ORF Size: 876 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026490.2](#), [NP_080766.1](#)

RefSeq Size: 5016 bp


RefSeq ORF: 879 bp

Locus ID: 56284

UniProt ID: [Q9D338](#)

Cytogenetics: 6 35.81 cM

Product images:



Circular map for MG204003