

Product datasheet for **MG204155**

Jam2 (NM_023844) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Jam2 (NM_023844) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Jam2
Synonyms:	1110002N23Rik; 2410030G21Rik; 2410167M24Rik; JAM-2; JAM-B; Jcam2; VE-JAM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204155 representing NM_023844 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGAGGAGCCCCAAGGCCTCCTGATGCTGCTGCTACACTACTTGATCGTCGCCCTGGACTATC
ATAAGGCAAATGGGTTTTCTGCATCAAAGACCACCGTCAAGAAGTCACAGTAATAGAGTTCCAAGAGGC
TATTTTGGCTTGTAAAACCCAAAGAAGACTACCTCCTCCAGACTGGAGTGAAGAAGGTGGACAGGGG
GTCTCCTTGGTCTACTACCAACAGGCTCTCCAAGGTGACTTTAAAGACCGTGTGAGATGATAGATTTCA
ATATACGAATCAAAAATGTTACAAGAAGTGATGCTGGAGAGTATCGCTGTGAAGTCAGCGCTCCGACTGA
GCAAGGCCAGAACCTGCAGGAAGATAAAGTCATGCTAGAAGTACTAGTGGCTCCTGCTGTTCTGCCTGT
GAAGTGCCCACTTCTGTTATGACTGGAAGTGTGGTGGAGTACGATGCCAGGATAAAGAAGGAAACCCAG
CTCCGGAGTACATCTGGTTAAAGATGGCACAAGTTTGTAGGGAATCCAAAAGGCGGCACACACAACAA
CAGCTCGTACACAATGAACACGAAGTCTGGAATTCTGCAATTCACATGATTTCCAAGATGGACAGTGGAG
GAGTATTACTGCGAAGCCCGAACTCTGTCGGACACCGCAGGTGCCCTGGGAAGCGAATGCAAGTAGATG
TTCTCAACATAAGCGGCATCATAGCAACGGTTGTGGTGGCCCTTCGTATTTCTGTATGTGGCCTTGG
CACATGCTATGCTCAGAGGAAAGGCTACTTTTCAAAGAACTTCCTTCCAGAAGGGCAGTCTGCATCT
AAAGTCACTACGATGAGCGAAAAATGATTTCAAGCACACAAAATCCTTTATAATT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG204155 representing NM_023844
 Red=Cloning site Green=Tags(s)


MARSPQGLLMLLLLHYLIVALDYHKANGFSASKDHRQEVTVIEFQEAILACKTPKKTSSRLEWKKVGQG
 VSLVYYQALQGDGFKDRAEMIDFNIRIKNVTRSDAGEYRCEVSAPTEQGQNLQEDKVMLEVLVAVPAC
 EVPTSVMTGGSVELRCQDKEGNPAPEYIWFKDGTSLLGNPKGGTHNNSYTMNTKSGILQFNMISKMDSG
 EYYCEARNVSVGHRRCPGKRMQVDVLNISGIIATVVVAVFVISVCGLGTCTYAQRKGYFSKETSFKGSPAS
 KVTMTSENFKHTKSFII

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja1824_g10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_023844


ORF Size: 894 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_023844.5 , NP_076333.3
RefSeq Size:	2902 bp
RefSeq ORF:	897 bp
Locus ID:	67374
UniProt ID:	Q9J159
Cytogenetics:	16 C3.3
Gene Summary:	May play a role in the processes of lymphocyte homing to secondary lymphoid organs. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG204155