

Product datasheet for **MG204290**

Lysmd3 (NM_030257) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	TurboGFP
Symbol:	Lysmd3
Synonyms:	1110030H10Rik; BC003322
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)

ORF Nucleotide Sequence: >MG204290 representing NM_030257
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCGGGGAGGAATCAGAACCGGACGGTTTCCCTCCCAGGAATTCAGGCCAGTGCCATGTACTTGCTT
TTGAAATTGTACAGACAATGACATGTTGGAAGAAGATGCTGAAGTCTACGAATTCGGTCCAGGGGAAA
GGAGAAGGTTCAAGAAGCGCATCGAGAGATAGACTTGATGACATTGTTATATTAACCAAAGATATCCAG
GAAGGAGACACTCTAAACGCAGTAGCCCTTCAGTACTGCTGTACGGTAGCAGATATCAAAAGAGTTAACA
ATCTCATCAGTGATCAAGACTTTTTGCTCTTAGGTCTATCAAAATTCAGTTAAAAGGTTTAGTTCCTTT
GACTGAACTCTTCATCCTCTGAAAGGAAGACACATTTACATCCTCCACCTGTTCCGTATTTTCAAGAG
CAAGACATTGTGCCAGCTGATGGTTCTTTCTCCAGTGAGTCAGCCGGCAGTTTCCTAAAAGAAGTGG
ACCGAGACATAGAACAATAGTCAAGGTACAGACACCAAGAAGGAGAACCTGAATGAAGTCGTGTCTGC
CTTAACAGCACAACAGGTCCGATTTGAACCTGATAACAAAAGCATTACCGAAAAGGATCCGTATTACGGA
GCAGACTGGGGATTGGTTGGTGGACAGCTGTAGTGATAATGTTGATAGTGGGTATAATAACACCAAGTGT
TTTATCTGCTGATTATGAAATTTAGCTAAAGTGGATGTTAGTCACCATTCACAGTGGGTTCTTCACA
CTTGATCCAGGACTCACACCTCCAACACAGCAGAGAGATGGAAATGAAATTTGGTCCAACAAAAGGA
ATACCTGTTGGTCAGCAAGATGACCACAACTATATAGGCAAGACCCTCAGGCACATGATGCTCAACACA
AAACG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



Protein Sequence: >MG204290 representing NM_030257
 Red=Cloning site Green=Tags(s)



MAGRNNQRTVSLPGIQASGHVLAFGNCTDNDMLEEDAIEVYELRSRGKEKVRSSASDRDLDDIVILTKDIQ
 EGDTLNAVALQYCCTVADIKRVNNLISDQDFALRSIKIPVKRFSSLTETLHPLKGRHILHPPVPYFQE
 QDIVPADGSLSSSESAGSFLKEVDRDIEQIVKCTDTKKENLNEVVSALTAQQVRFEPDNKSIHRKDPYYG
 ADWGIGWWTAVVIMLIVGIITPVFYLLYYEILAKVDVSHHSTVGSSHLHPGLTPPTQHREMENEIGPTKG
 IPVGQQDDHKLYRQDPQAHAQHKT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_030257

ORF Size: 915 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_030257.1](#), [NP_084533.1](#)

RefSeq Size: 3139 bp

RefSeq ORF: 918 bp


Locus ID: 80289

UniProt ID: [Q99LE3](#)

Cytogenetics: 13 C3

Gene Summary: Essential for Golgi structural integrity.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG204290