

## Product datasheet for **MG204297**

### 1810043G02Rik (NM\_026431) Mouse Tagged ORF Clone

#### Product data:

|                           |                         |
|---------------------------|-------------------------|
| Product Type:             | Expression Plasmids     |
| Tag:                      | TurboGFP                |
| Symbol:                   | 1810043G02Rik           |
| Synonyms:                 | AV026620; D10Jhu13e     |
| Mammalian Cell Selection: | Neomycin                |
| Vector:                   | pCMV6-AC-GFP (PS100010) |
| E. coli Selection:        | Ampicillin (100 ug/mL)  |

**ORF Nucleotide Sequence:** >MG204297 representing NM\_026431  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGAAGCTGACACGAAAGATGGTCCTGTCCCGAGCCAAGGCTTCGGAGCTGCACAACGTGCGGAAGCTCA  
ACTGCTGGGGCAGCCAGCTCACAGATATTTCCATATGCCGGGAGATGCCCAGCCTAGAGGTAATCACACT  
CAGTGTCAACAGTGTCTCCACCCTGGAACCGGTGCGCAGCTGCAGGCGCCTGAGCGAGCTGTACCTGAGG  
AGGAACCGCATCCCGAGTCTGAACGAGCTGTTCTACCTGAAAGACCTGCCACACCTCAGGGTGCTCTGGC  
TGGCCGAGAACCCGTGCTGTGGCACCAGCCCGCACCTTTACCGCATGACCGTGCTGCGCAACCTGCCCA  
CCTGCAGAAGCTAGACAACCAAGCTGTGACTGAGGAAGAGCTGACCCGCGCGCTGATGGAGGGAGATGAG  
ATCACTGCCGCCCCGCACAGAGAAGGTGCAGGCAATGGCTGTCCCAAGCCCCCTACGCCCTCCACTCTG  
TCAGCTCTGCAACTGAGACCAGCCAGCACCTGCTGAGCTACACAGAGGAGACAGAAGTCCAGGGCCAGAC  
TACCACAGACCAGTCTCCTTCTTCTACCAAGGGACACGATGAGGAGTCACAAGAGTAGGAACATCCTG  
ACTGCCATATTGCTACTGCTGCGGGAGCTGGATACGGAGGGGCTGGAGACTGTCCAGCAGACCGTGGGCA  
GCCGGCTACAGGCGCTGCACAGGCCGAGCCCCAGGAAGACATGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



**Protein Sequence:**

>MG204297 representing NM\_026431

Red=Cloning site Green=Tags(s)


MKLTRKMLVSRKASELHNVRLNCWGSQSLTDISICREMPSEVITLSVNSVSTLEPVRSCRRLESELYLR  
RNRIPSLNELFYKDLPHLRVLWLAENPCCGTSPHLYRMTVLRNLPHLQKLDNQAVTEEEELTRALMEGDE  
ITAAPHREGAGNGCPKPPYALHSVSSATETSQHLLSYTEETEVEVGGTTTDDQSPSPFSPRDTMRSHKSRNII  
TAILLLLRELDTEGLETVQQTVGSRLQALHRPEPQEDME

TRTRPLE - GFP Tag - V

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**ACCN:**

NM\_026431

**ORF Size:**

747 bp

**OTI Disclaimer:**

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:**


This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:**

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <u>NM_026431.2, NP_080707.2</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq Size:</b>           | 1696 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 750 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 67884                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UniProt ID:</b>            | <u>Q8C6G1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 10 39.72 cM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Gene Summary:</b>          | Plays a role in cilia formation and/or maintenance (PubMed:21289087). Plays a role in the regulation of cell morphology and cytoskeletal organization (By similarity). Involved in DNA damage repair (By similarity).[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                                                                                                    |

### Product images:



Circular map for MG204297