

Product datasheet for **MG204314**

Ung (BC011039) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ung (BC011039) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ung
Synonyms:	UNG1, UNG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204314 representing BC011039 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATCGGCCAGAAGACCCTATACTCTTTCTTCTCTCCGACCCACAGGGAAGCGGACCACGCGCAGCC
CGGAGCCGGTCCCAGGGAGCGGAGTGCCGGCCGAGATCGGCGGCGATGCCGTGGCCAGCCCCGCAAGAA
GGCCCCGGTGGAGCAGAACGAGCAAGGCTCGCCGCTCAGCGCCGAGCAGCTCGTCCGCATCCAGAGGAAC
AAGGCCGCGCGCTGCTCAGGCTCGCCGCCGCAACGTGCCTGCCGGCTTCGGCGAGAGCTGGAAGCAGC
AGCTGTGCGGGGAGTTCGGGAAGCCGTAATTCGTCAAGCTAATGGGATTTGTTGCTGAAGAAAGGAACCA
CCACAAGGTCTATCCGCCCCGGAGCAGGTGTTACATGGACCCAGATGTGCGACATCCGAGATGTGAAG
GTTGTCATTCTGGGACAGGATCCCTATCACGGACCTAATCAAGCTCACGGGCTCTGCTTCAGTGTCCAAA
GACCAGTTCACCCCCGCCAGTTTGGAAAACATTTTTAAAGAGCTGTCTACAGACATCGATGGTTTTGT
TCATCCTGGCCATGGGGATTTGTCAGGGTGGGCCCGACAAGGTGTCTCCTCCTCAACGCCGCTCCTCACT
GTCCGCGCCACCAAGCCAATCCCATAAGGAGAGGGGCTGGGAGCAGTTCACGGACCGGTCGTGTCTCT
GGCTGAATCAGAACCTGAGTGGCCTCGTCTTCTCCTCTGGGGCTCTTATGCTCAGAAGAAGGGCAGTGT
CATTGACAGGAAGCGTCACCATGTTCTGCAGACAGCTCACCCCTCCCCGTGTCGGTGCACAGAGGGTTC
CTGGGATGCAGACATTTTTCTAAAGCCAATGAAGTGTCCAGAAGTCTGGCAAGAAGCCCATCAACTGGA
AGGAGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG204314 representing BC011039
 Red=Cloning site Green=Tags(s)

MIGQKTLYSFFSPTPTGKRTRRSPEVPVPGSGVAEIGGDAVASPAKKARVEQNEQGSPLSAEQLVRIQRN
 KAAALLRLAARNVPAGFGESWKQLCGEFGKPYFVKLMGFVAERNHHKVYPPEQVFTWTQMCDIRDVK
 VVILGQDPYHGNQAHGLCFVQRPVPPPPSLENIFKELSTDIDGFVHPGHGDLSGWARQGVLLNNAVLT
 VRAHQANSHKERGWEQFTDAVSWLNQNL SGLVFLWGSYAQKKGSVIDRKRHHVLQTAHPSPLSVHRGF
 LGCRRHFSKANELLQKSGKKPINWKEL


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: BC011039

ORF Size: 920 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC011039 , AAH11039
RefSeq Size:	2002 bp
RefSeq ORF:	920 bp
Locus ID:	22256
Cytogenetics:	5 F
Gene Summary:	Excises uracil residues from the DNA which can arise as a result of misincorporation of dUMP residues by DNA polymerase or due to deamination of cytosine.[UniProtKB/Swiss-Prot Function]