

Product datasheet for **MG204414**

Krt13 (BC060286) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Krt13 (BC060286) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Krt13
Synonyms:	CK13, K13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204414 representing BC060286 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCTGCCGCTTCCAGAGTTCCTCCATGAGCTATGGAGGTGGTTTCGGAGCTGGTTCTTGCCAGCTTG
GAGGAGGCCGTAATATCTCCTCCTGCTCATCTCGGTTTGTCACTGGAGGGTCAGCTGGAGGCTATGGAGG
GGGCATGAGCTGTGGCTTTGGTGGAGGAGCTGGTGGTGGCTTTGGAGGTGGCTTTGGAGGTGGCTTTGGA
GGTAGCTATGGCGGAGGTTTGGCGGTGGCTTTGGTACTTCGGAGGTGTCGATGGTGGCCTCCTCTCTG
GCAACGAGAAGATCACCATGCAGAACCTCAATGACCGCCTGGCCAGCTACCTGGATAAAGTGGCGCCTT
GGAGGCAGCCAATGCCGACCTGGAGGTGAAGATTCGTGACTGGCATCTGAAACAGAGCCAGCTAGTCCA
GAGCGGGACTACAGTGCTTATTATAAGACCATCGAGGAGCTCCGGATCAAGATTCTGGAAGCCACCACGG
ACAACAACCGGATCATCCTGGAGATTGACAATGCCAGGCTGGCTGCAGATGACTTCAGGCTCAAGTACGA
AAATGAGCTGACCTTGCGCCAGAGCGTGGAGGCCGACATCAACGGCCTGCCAGGGTGTGGATGAGCTG
ACTCTGGCTAAGACTGACCTGAAATGCAGATTGAGAGCCTGAATGAAGAAGTGGCTTACCTGAAGAAGA
ACCACGAAGAGGAGATGAAGGAATTCAGCAACCAGGTAGTAGGCCAGGTCAATGTGGAGATGGACGCCAC
CCCTGGCATTGATCTGACCCGTGTGCTGGCTGAGATGAGGGAGCAGTATGAAGCCCTGCCAGAGAAGAAT
CGGAGGGATGCTGAGGAATGGTTCCAGACCAAGAGTGCAGAGCTGAACAAGGAAGTATCCTCCAATTCCG
GTCGCCAGATTTCCGAAAGTAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG204414 representing BC060286
 Red=Cloning site Green=Tags(s)

```


MSCRFSQSSMSYGGGFGAGSCQLGGGRNISSCSSRFVTGGSAGGYGGGMSCGFGGGAGGGFGGGFGGGF
GSYGGGFGGGFDGFGVDGGLLSGNEKITMQNLNDRLASYLDKVRALEAANADLEVKIRDWHLKQSPASP
ERDYSAYYKTIEELRIKILEATTDNNRIILEIDNARLAADDFRLKYENELTLRQSV EADINGLRRVLDL
TLAKTDLEMQIESLNEELAYLKKNH EEMKEFSNQVVGQVNVEMDATPGIDLTRVLAEMREQYEALAEKN
RRDAEEWFQTKSAELNKEYSSNSGRPDFRKY
  
```

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: BC060286

ORF Size: 933 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC060286.1
RefSeq Size:	1232 bp
RefSeq ORF:	935 bp
Locus ID:	16663
Cytogenetics:	11 63.41 cM
Gene Summary:	The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. This type I cytokeratin is paired with keratin 4 and expressed in the suprabasal layers of non-cornified stratified epithelia. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]