

Product datasheet for **MG204676**

Akr1c14 (NM_134072) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Akr1c14 (NM_134072) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Akr1c14
Synonyms: 9030611N15Rik; AW557061
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG204676 representing NM_134072
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATTCTGTATCCCCACGTGTGGTACTAAACGATGGTCACTTCATCCCTGCACTGGGGTTTGGAACCA
CTGTGCCCGATAAGGTTCTAAGGATGAAGTATCAAGGCTACTAAAATAGCTATAGATACTGGATTCCG
CCATTTTGACTCCGTTATTTGTACCAATAGAAGAGGAAGTAGGCCAGGCCATTAGAAGCAAGATTGAA
GATGGCACTGTGAAGAGAGAAGATATTTCTATACTTCAAAGCTTTGGAGCACTTCCATCGACCAGAAT
TGGTCCGATCTTGCTTGGAAAAGACACTGAAGAATGCACAACCTGGACTATGTTGATCTTTACATTATTCA
TTTCCCAATGGCTTTGCAGCCTGGAGATAAACTATTTCCACGAGACGAACATGGAAAACCTTTGGCTGAA
GCAGTGGATCTCTGTGACACATGGGAGGCCATGGAAAAGTGTAAAGATGCCGGACTGGCCAAGTCCATCG
GAGTGTGCAACTTTAACTTCAGGCAGCTGGAGACGATTCTGAACAAGCCGGGGCTCAAGTACAAGCCTGT
GTGCAACCAGGTAGAATGCCATCTTTATTTAAACCAGAGCCAAATGCTGGACTATTGTAAGTCAAAAAGAC
ATCATTCTGGTTTCTACTGCACATTGGGAAGTTCACGAGACAAAATCTGGGTGGACCAGAAAAGTCCAG
TTCTCTTAGATGATCCAGTTCTTTGTGCCATGGCAAATAAGTACAAGCAAACACCAGCACTGATTGCCAT
TCGTTACCAATTACAGCGTGAATTGTGGTCTGACCAGGAGTTTCAAGGAGAAGCGGATCAAAGAGTTC
ATGAAGTTTTTGAATTCCAGTTGGCTTCAGAGGACATGAAAGTCTGGATGGCTTGCACAGAAATTTAA
GATACAATACTGCGAGTTATTTTGATGACCATCCCAATCATCCATTTACTGATGAATAT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG204676 representing NM_134072
Red=Cloning site Green=Tags(s)

MNSVSPRVVLDNGHFIPALGFGTTVPDKVPKDELIKATKIAIDTGRHFDSAYLYQIEEEVGVQAIRSKIE
 DGTVKREDIFYTSKLWSTFHRPELVRSCKLEKTLKNAQLDYVDLYIIHFPMALQPGDKLFPRDEHGKLLAE
 AVDLCDTWEAMEKCKDAGLAKSIGVSNFNFRQLETILNKPGLKYKPVCNQVECHLYLNQSQMLDYCKSKD
 IILVSYCTLGSSRDKIWVDQKSPVLLDDPVLCAMANKYKQTPALIAIRYQLQRGIVVLRFSFKEKRIKEF
 MKVFEFQLASEDMKVL DGLHRNRLRYNTASYFDDHPNHPFTDEY

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_134072

ORF Size: 969 bp

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_134072.2 |
| RefSeq Size: | 2459 bp |
| RefSeq ORF: | 972 bp |
| Locus ID: | 105387 |
| Cytogenetics: | 13 A1 |