

Product datasheet for **MG204726**

Ccdc49 (BC079553) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ccdc49 (BC079553) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ccdc49
Synonyms:	1300013D05Rik; AA473636; Ccdc49; R75228
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204726 representing BC079553 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAGGCGGAGACCTGAACCTGAAGAAGAGCTGGCACCCGAGACCCTTCGCAATGTGGAGAAAGTGT
GGAAGGCGGAGCAGAAGCATGAAGCGGAGAGGAAGAAGATTGAGGAGTTCAGCGAGAGCTCCGAGAGGA
GAGAGCGGGGAGGAGATGCAGCGCTACGCTGAGGATGTTGGGGCTGTCAAGAAAAAGAAGAAAGTTG
GATTGGATGTACCAAGGTCGGGTGGCATGGTGAACCGTGATGAGTACTTGTGGGGCGCCCCATTGACA
AGTATGTTTTTGAAGAATGGAGGAGAGGAAGCTGGCTGCTTCTGAGACAGGACTCCTCCCAGGCTC
AATCTTTGCCCGTCGGGCGCAATTCCTTCTTGACATGGCCAGCAAAATCCGGGAAGACCTCTCTTC
ATCATAAGGAAGAAAGAGGAGGAGAAAAAAGAGAAGTGCTGAACAATCCTGTGAAAATGAAGAAAATCA
AAGAGCTGTTGCAATGAGCCTTGAGAAAAAGGAGAAGAAGAAGAAGGAGAAGAAGAAGACAG
GAAGCACAAACACAGAAGCTCCAGCAGCGGTGGCTCCAGCAGCGAGGATGAGCAGAGCCAGGCCCGGTCT
CAAAGAAGATGGCAAATTCCTTCTGTTTTGTCAAAGTCCCTGGATATGGCTTACAGGTGAGGGACT
CTGACCGAACCAGGGGCTGCAGGTTCTTGGGTGAACAGAGGGCGATAAAGAACAATCCCGGTCAAG
AAGCTCCTCACCCCCAGGCATGCCAGCAAGAAGAGCACCAAGGAAGAGAGGCCCGGGACAGGAGGTCT
CGATCCCCCAGCAGAAGTCCCGGTCCCAAGGCCAGCAAACCGCATACCTCTAAGGTAACAGGAAAG
AGAGAGACAGCCCATCACCTAAGAAGGAGGCCTACCAGAGCGGCATGCTTCTGGGTACACCAGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG204726 representing BC079553
 Red=Cloning site Green=Tags(s)

MGGDDLNLKKSWHPQTLRNVEKVKAEQKHEAERKKIEELQRELREERAREEMQRYAEDVGAVKKKEEKL
 DWMYQGP GGMVNRDEYLLGRPIDKYVFEKMEEREAGCSSETGLLPGSIFAPSGANSLDMASKIREDPLF
 IIRKKEEEKKREVLNNPVKMKIKELLQMSLEKKEKKKKKKKKHRKHKHRSSSSGGSSSEDEQSQARS
 QKKMANSFPVLSKVPGYGLQVRSDRNRLQGSLEQRAIKNNSRSRSPRPHASKKSTKEERPRDRRS
 RSPSRRSRSPRPSKPHTSKVNRKERDPSPKKEAYQRRHASGYTR

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC079553

ORF Size: 977 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC079553](#), [AAH79553](#)

RefSeq Size: 5862 bp


RefSeq ORF: 977 bp

Locus ID: 67480

Cytogenetics: 11 D

Gene Summary: Involved in pre-mRNA splicing as component of the spliceosome.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG204726