

Product datasheet for **MG204792**

4930430F08Rik (NM_175128) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	4930430F08Rik (NM_175128) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	4930430F08Rik
Synonyms:	4930571E09Rik; A730088G13Rik; AW551239; C430008C19Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204792 representing NM_175128 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCGCTGGGCTCGGTGCAGCGGAAGATGCCGTGTGTGTTGTGACGGAGGTGAAAGCTGAGCCTT
CAGCCAAGCGGGAGCATCAGCCATTTAAAGTTTTGGCAACTGAACTCTAAGTGAAAAGGCGTTAGATGC
AGACGTGTACAATGCGGTTGCAACAGAAAAAGTGGATGGAACATGCTGCTATGTTACTAACTACAAAGT
CAGCCATACCTTTGGGCTCGACTAGATAGAAAACCTAACAAAGCAGGCTGACAAGAGATTTAAAAATTTT
TACATTTAAAAGAAAGTGCAAAAGAGTTTCACTGGAACACTGAGGAAGACTTCAAGCCTGTCCAGAGTG
CTGGATACCAGCAAAGGAAATAGAAAAACAAAACGGGAAGCCAGTGCCTGATGAAAACGGACACATCCCT
GGTTGGGTACCAGTAGAGAAGGGCAGCAAGCAATATTGCTGGCACTCGTCGGTAGTGAATTACGAGTTTG
GAATTGCCCTGGTACTAAGGCACCATCCTGATGATCCCGGGTTTTGGAAATTAGTGCTGTTCCCTGTG
AGAGCTTTTAGAGCAAACCCTGGAGCTGATAGGCACGAGTATCAATGAAAACCCCTATGGGCTGGGAAGC
AAGAAATCCCGTTACATTTTCTCACACCACACGGAGCATTCCAGGTGAGAAATCTGCCACACTGAAAC
ACAATGACCTGCTGTCTGGTTGAAGACTGCAGAGAGGGTCAAATTGAAGGGATTGTGTGCATTGCGG
TGACGGCTGTTTAATTAAGGTCCATCGTCACCATCTTGGTTTATGCTGGCCACTCCAGATACTTACATG
AATTCAAACCAGTTATTATTAACATGAACCTGAACCTGAACAATTATGACTGTCTTTTGATAACCACT
CTTTGTTAATCAGTTTTCAAAAATAGATAAACAGAAATTTGAGAGACTCAAAGACATAAATACTTGATGT
A

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG204792 representing NM_175128
Red=Cloning site Green=Tags(s)

MKRLGSVQRKMPQVFVTEVKAEPSAKREHQPFKVLATETLSEKALDADVNAVATEKVDGTCCYVTNYKG
 QPYLWARLDRKPNKQADKRFKFLHSKESAKEFHWNTEEDFKPVPECWIPAKEIEKQNGKPVDPENGHIP
 GWVPVEKGSQYCWSSVYNEFGIALVLRHHPDDPGVLEISAVPLSELLEQTLELIGTSINGNPYGLGS
 KKSPLHFLTPHGAFQVRNLPTLKHNDLLSWFEDCREGQIEGIVWHCGDGLIKVHRHHLGLCWPLPDTYM
 NSKPVIINMNLNLYDCAFNDQSLFNQFSKIDKQKFERLKDIIIDV


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_175128

ORF Size: 981 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_175128.2
RefSeq Size:	2375 bp
RefSeq ORF:	984 bp
Locus ID:	68281
UniProt ID:	Q8BHN7
Cytogenetics:	10 D1