

Product datasheet for **MG204927**

Gapdh (BC091768) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gapdh (BC091768) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gapdh
Synonyms:	Gapd
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG204927 representing BC091768 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGAAGGTCGGTGTGAACGGATTTGGCCGTATTGGGCGCCTGGTCACCAGGGCTGCCATTTGCAGTG
GCAAAGTGGAGATTGTTGCCATCAACGACCCCTTCATTGACCTCAACTACATGGTCTACATGTTCCAGTA
TGACTCCACTCACGGCAAATTCACGGCACAGTCAAGGCCGAGAATGGGAAGCTTGTCATCAACGGGAAG
CCCATCACCATCTTCCAGGAGCGAGACCCCACTAACATCAAATGGGGTGAGGCCGGTCTGAGTATGTCCG
TGGAGTCTACTGGTGTCTTACCACCATGGAGAAGGCCGGGCCACTTGAAGGGTGGAGCCAAACGGGT
CATCATCTCCGCCCTTCTGCCGATGCCCCATGTTTGTGATGGGTGTGAACCACGAGAAATATGACAAC
TCACTCAAGATTGTCAGCAATGCATCCTGCACCACCAACTGCTTAGCCCCCTGGCCAAGGTCATCCATG
ACAACCTTTGGCATTGTGGAAGGGCTCATGACCACAGTCCATGCCATCACTGCCACCCAGAAGACTGTGGA
TGGCCCTCTGGAAAGCTGTGGCGTGATGGCCGTGGGGCTGCCCAGAACATCATCCCTGCATCCACTGGT
GCTGCCAAGGCTGTGGCAAGGTCATCCAGAGCTGAACGGGAAGCTCACTGGCATGGCCTTCCGTGTTCT
CTACCCCAATGTGTCCGTCTGGATCTGACGTGCCGCTGGAGAAACCTGCCAAGTATGATGACATCAA
GAAGGTGGTGAAGCAGGCATCTGAGGGCCCACTGAAGGGCATCTTGGGCTACACTGAGGACCAGGTTGTC
TCCTGCGACTTCAACAGCAACTCCCACTTCCACCTTCGATGCCGGGGCTGGCATTGCTCTCAATGACA
ACTTTGTCAAGCTCATTTCCTGGTATGACAATGAATACGGCTACAGCAACAGGGTGGTGGACCTCATGGC
CTACATGGCCTCAAGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG204927 representing BC091768
 Red=Cloning site Green=Tags(s)

MVKVGVNGFGRIGRLVTRAAICSGKVEIVAINDPFIDLNYMVYMFQYDSTHGKFNGTVKAENGLVINGK
 PITIFQERDPTNIKWGEAGA EYVVESTGVFTTMEKAGAH LKGGAKRVII SAPSADAPMFVMGVNHEKYDN
 SLKIVSNASCTTNCLAPLAKVIHDNFGIVEGLMTTVHAI TATQKTVDGPSGKLWRDGRGAAQNIIPASTG
 AAKAVGKVIPELNGKLTGM AFRVPTPNVSVVDLTCRLEKPAKYDDIKKVVKQASEGPLKGI LGYTEDQVV
 SCDFNSNSHSSTFDAGAGIALNDNFVKLISWYDNEYGYSNRVV DLMAYMASKE


TRTRPLE - GFP Tag - V

Restriction Sites:


Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: BC091768

ORF Size: 999 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC091768.1
RefSeq Size:	1237 bp
RefSeq ORF:	1001 bp
Locus ID:	14433
Cytogenetics:	6 59.32 cM
Gene Summary:	This gene encodes a member of the glyceraldehyde-3-phosphate dehydrogenase protein family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. The encoded protein was originally identified as a key glycolytic enzyme that converts D-glyceraldehyde 3-phosphate (G3P) into 3-phospho-D-glyceroyl phosphate. Subsequent studies have assigned a variety of additional functions to the protein including nitrosylation of nuclear proteins, the regulation of mRNA stability, and acting as a transferrin receptor on the cell surface of macrophage. Alternative splicing results in multiple transcript variants. Many pseudogenes similar to this locus are found throughout the mouse genome. [provided by RefSeq, Jan 2014]