

## Product datasheet for **MG205009**

### Zscan22 (BC026676) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Zscan22 (BC026676) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Zscan22
Synonyms:	AW227621; D530006B18Rik; Hkr2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205009 representing BC026676 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGCGACTGAACTGTCATGGGAAGCGTTTTGCCAAAATCTACCCTTGCCACACTTGTAACCTGAGA  
GCCACTCAGAGAGCCAGCCAGAACTCCTAGGAGCACTCTGGATGAAGTCTACTGCCAAGAGATGGATTT  
CAGAAAAGCTCTGGGGCCTCACATGGATGCCCCAAAGACCAGCCTGGCCATGAGTCTAACACCTCAGGA  
AACGGTCCAACATGTGGCCCACTTCCCTTCCAAGACAAGGCTTCTTCAGAGGAGAAAATTTGGCCCAT  
TACTTGATAATGAGACAGTGCCTCCAGATACTTGCTCAGAGAAGAAATCCTCCAAGGATAGTGAATGTTT  
GAAAACATTCCAGAACACTTCAGCTCTGGAGGCCACCAGAAGAGCCACTCTCAGAAGACGCCCTATGCC  
TGTAAGGACTGCGTGTAAAGAGTGTGGTAAAGCCTTCAGTCGCGTTGCCACCTGACACAGCACCTGAG  
GATTCATACTGGGGAGAAGCCTTATAAATGTGAGGAATGTGGTAAAACCTTCAGCCGAGCACCCACCTC  
ACGCAGCACCCAGCGAGTACACACAGGCGAGCGCCCTACGAGTGTGATACTTGCGGCAAAGCCTTCAGCC  
AGAGCACCCACCTGACACAGCACCCAGCGCATCCACACCGGGGAGAAGCCCTACAGGTGTGATGTGTGG  
GAAAGCCTTCAGTGACTGTTTCAGCTCTGGTCCGGCACCTGAGAGTCCACTCTGGAGAGAAGCCTTATCAA  
TGTAAGGACTGCCGAAGCCTTCGCACAAAGCTCCTCCCTCATTGAGCACCCAGCGCACTCACACAGGTG  
AAAAGCCTTATAAATGCAGTGACTGTGGGAAAGCTTTTAGCCGCAGCTCAGCTCTCATGGTGCACCTAAA  
GATCCATATCACCGTGACTCGAAGCACCCCT

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG205009 representing BC026676  
 Red=Cloning site Green=Tags(s)

MATETVMGSVLPKSTLAHTCKPESHSESQPELLGALWMKSTAQEMDFRKALGPHMDAPKDQPGHESNTSG  
 NGSNMWPNFPSQDKASSEKFGPLLDNETVPPDTCSEKSSKDSECLKTFQNTSALEAHQKSHSQKTPYA  
 CTECGKVF SRSTHLVQHQQVHTGAKPHACKECGKAF SRVAHLTQHRLRIHTGEKPYKCEECGKTF SRSTHL  
 TQHQ RVHTGERPYECDTCGKAF SQSTHLTQHQR IHTGEKPYRCDVCGKAFSDCSALVRHLRVHSGEKPYQ  
 CKDCPKAFAQSSSLIEHQRTHTGEKPYKSDCGKAF SRSSALMVHLKIHTVTRSTP

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** BC026676

**ORF Size:** 1013 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC026676</a> , <a href="#">AAH26676</a>
<b>RefSeq Size:</b>	2124 bp
<b>RefSeq ORF:</b>	1013 bp
<b>Locus ID:</b>	232878
<b>Cytogenetics:</b>	7 A1