

Product datasheet for **MG205041**

Fnbp1 (NM_019406) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fnbp1 (NM_019406) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Fnbp1
Synonyms: 1110057E06Rik; 2210010H06Rik; FBP1; Fbp17
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG205041 representing NM_019406
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAGCTGGGGCACTGAGCTCTGGGATCAGTTTGACAACCTGGAAAAACATACACAGTGGGGAATCGATA
TTCTTGAGAAATACATCAAGTTTGTCAAGGAGAGGACGGAGATTGAGCTCAGCTATGCCAAGCAACTCAG
GAATCTTTCAAAGAAATACCAACCTAAGAAGAACTCGAAGGAAGAGGAGGAGTACAAGTACACGGCTTGC
AAAGCCTTTCTTCCACCCTGAATGAGATGAATGACTACGCCGGGCAGCAGGAGTATCTCTGAGAACA
TGACGTCACAGATCACGGTGGACCTGATGCGCTACGTTTCAGGAGCTGAAGCAGGAGAGGAAATCGAACTT
CCATGATGGACGGAAGGCTCAGCAGCACATAGAAACGTGTTGGAAGCAACTGGAGTCAAGTAAGAGGAGG
TTTGAGCGGGACTGTAAGGAAGCCGACCGGGCACAGCAGTACTTCGAGAAAATGGACGCTGACATCAACG
TGACCAAGGCGGATGTGGAAAAGGCACGACAACAAGCTCAGATACGCCAGCAAAATGGCAGAGGACAGCAA
AGCAGATTACTCCTTAATCCTGCAGAGATTCAACCAGGAGCAATGGGAATACTACCATACCCACATCCCC
AACATCTTCCAGAAAATACAAGAGATGGAGGAGAGGGGATTGTGCGTATTGGGGAGTCCATGAAGACGT
ACGCAGAGGTGGACCGGCAGGTGATACCCATCATCGGGAAGTGCCTGGACGGGATAGTGAAGCGGCCGA
GTCTATCGACCAGAAAAACGACTCCCAGCTGGTTCGTAAGCCTATAAGTCAGGATTCGAGCCTCCTGGA
GACATTGAATTCGAAGATTACACACAGCCAATGAAACGCACAGTGTGACACAACAGCCTTCCAGCTCCA
AAGAAGGCAAGCCTGAGCTCAGATTTGGCGGCAAGTCCAGAGGCAAGCTCTGGCCATTATCAAGAAAAA
CAAGTACTGGCCATTTGGACCCTGCGTGGGCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG205041 representing NM_019406
Red=Cloning site Green=Tags(s)

MSWGTELWDQFDNLEKHTQWGIDILEKYIKFVKERTEIELSYAKQLRNL SKKYQPKKNSKEEEEEYKYTAC
 KAF LSTL NEMNDYAGQHEVISENMTSQITVDLMRYVQELKQERKSNFHDGRKAQQHIETCWKQLESKR
 FERDCKEADRAQQYFEKMDADINVTKADVEKARQQAQIRQMAEDSKADYSLILQRFNQEQWEYYHTHIP
 NIFQKIQEMEERRIVRIGESMKTYAEVDRQVIPIIGKCLDGI VKAAESIDQKNDSQLVVEAYKSGFPEPPG
 DIEFEDYTQPMKRTVSDNSLSSSKEGKPELRFGGKSRGKLWPF IKKNKVLAIWTLRGL

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_019406

ORF Size: 1014 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019406.3 , NP_062279.1
RefSeq Size:	1977 bp
RefSeq ORF:	1017 bp
Locus ID:	14269
UniProt ID:	Q80TY0
Cytogenetics:	2 B
Gene Summary:	Required to coordinate membrane tubulation with reorganization of the actin cytoskeleton during the late stage of clathrin-mediated endocytosis. Binds to lipids such as phosphatidylinositol 4,5-bisphosphate and phosphatidylserine and promotes membrane invagination and the formation of tubules. Also enhances actin polymerization via the recruitment of WASL/N-WASP, which in turn activates the Arp2/3 complex. Actin polymerization may promote the fission of membrane tubules to form endocytic vesicles. May act as a link between RND2 signaling and regulation of the actin cytoskeleton. May be required for the lysosomal retention of FASLG/FASL (By similarity).[UniProtKB/Swiss-Prot Function]