

Product datasheet for **MG205247**

Cd72 (BC003824) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cd72 (BC003824) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cd72
Synonyms:	CD72c; Ly-19; Ly-32; Ly-m19; Lyb-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205247 representing BC003824 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGACGCTATCACGTATGCAGACCTGCGCTTTGTGAAAGTGCCCTGAAGAACAGTGCATCTAACC
ATCTAGGACAGGACTGTGAGGCCTATGAAGATGGGAACTCACCTACGAGAACGTGCAAGTGTCTCCAGT
CCCAGGAGGGCCACCAGGCTTGGCTTCCCCTGCACTAGCGGACAAAGCAGGGTTCGGGTGAGCAACCA
ACTGCGACCTGGAGCTCTGTGAAGTCGTCTGCTCTCAGGCAGATTCCCCGCTGCTCCTACAGTCTGCTTGC
AATACTTCTTGGCTTGGCCTTCTCGTGTCTGTCTGATGTTAGGGTGGCTGTCATCTGCCTGGGAGTTTCG
CTATCTGCAGGTGCTCGGCAGTTCAGGAGGGGACCAGGATTTGGGAAGCCACCAATAGCAGCCTGCAG
CAGCAGCTCAGGGAGAAGATAAGTCAGCTGGGGCAGAAGGAGGTGGAGCTTCAGAAGGCTCGGAAAGAGC
TGATCTCGAGCCAGGACACATTACAGGAGAAGCAGAGGACTCACGAGGACGCTGAGCAGCAACTACAAGC
CTGCCAGGCTGAGAGAGCGAAGACCAAGGAGAACCTGAAAACCTGAGGAGGAGCGGAGGAGGGACCTGGAC
CAGAGTTGACAAGCACGCGGGAGACACTAAGGCGCTTCTTCTGATTCATCAGACACCTGCTGTCCAT
GCGGATGGATTCCATATCAGAAAAGGTGCTTTACATCTCACATAACCTCAGAAGTCTGGAGGAAAGCCA
AAAATACTGCACATCTCTGTCTCCAACTGGCAGCATTCAATGACTCTCCTAAGTATTCCTATCAAGTT
TCTCTGCCAGCGGCTTAGAGGAGTTGCTAGATCGTTCGAAGTCATATTGGATACAGGAGATGAGCAAGA
AGTGGAGCTCATATTGTGACATGATACAAATATAACAAGAAGTGGGAAAGAACAATTTCCAAGTGTGC
AGAGCTTCACCCCTGCATTTGTGAGTTGGAGGCTTTCAGGTTTCTGATGGGATCCATCTGAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG205247 representing BC003824
 Red=Cloning site Green=Tags(s)

MADAITYADLRFVKVPLKNSASNHLGQDCEAYEDGELTYENVQVSPVPGPPGLASPALADKAGVGSEQP
 TATWSSVKSSALRQIPRCPTVCLQYFLLGLLVSLMLGVAVICLVGVRYLQVSRQFQEGTRIWEATNSSLQ
 QQLREKISQLGQKEVELQKARKELISSQDTLQEKQRTHEAEQQLQACQAERAKTKENLKTEEERRRDL
 QRLTSTRETLRRFFSDSSDTCCPCGWIPYQKRCFYISHNLRSLSEESQKYCTSLSSKLAAFNDSPKYSYQV
 SLPSPGLEELLDRSKSYWIQEMSKKWSSYCDMIQIYNKKWERTISKCAELHPCICELEAFRFPDGIHLN

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: BC003824

ORF Size: 1046 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC003824 , AAH03824
RefSeq Size:	1423 bp
RefSeq ORF:	1046 bp
Locus ID:	12517
Cytogenetics:	4 23.04 cM
Gene Summary:	Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]