

## Product datasheet for **MG205333**

### Ltb4r1 (NM\_008519) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ltb4r1 (NM\_008519) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ltb4r1  
**Synonyms:** BLT1; BLTR; Ltb4r; mBLTR  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG205333 representing NM\_008519  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGCTGCAAACACTACATCTCCTGCAGCACCTTCTTCTCCTGGTGGCATGTCCCTGTCTCTGTTGCCCA  
 TTGTTCTACTGTCTGTGGCCCTGGCTGTGGGGTCTCCTGGCAATAGCTTTGTAGTGTGGAGCATCTGAA  
 AAGGATGCAGAAACGCACGGTCACCGCCCTGCTGGTGTGAAGTGGCCCTGGCCGACTGGCTGTGTTG  
 CTCACTGCTCCCTTTTCTCCACTTTCTGGCTCGAGGCACCTGGAGTTTAGAGAGATGGGTTGCCGCC  
 TGTGCCACTATGTCTGTGAATAAGCATGTATGCCAGTGCCTGCTTATAACCATCATGAGTCTGGACCG  
 ATCACTGGCAGTGGCCCGCCCTTATGTCCAAAAGGTTGCTACTAAGGCCCTTGGCCGATGGGTGCTG  
 GCAGGCATCTGGGTGGTGTCTTTCTGCTGGCCATACCGGTCCTGTGTACCGTACAGTAAAATGGAACA  
 ACAGGACTCTGATCTGCGCTCCGAACTATCCCAACAAAGAGCATAAAGTCTTCCATCTGCTCTTCGAAGC  
 CATCACGGGCTTCTGCTGCCCTTCTAGCGGTGGTGGCCAGCTACTCTGACATCGGGCGCAGGCTGCAG  
 GCTCGGCGCTTCCGCCGTAGTCGCCGACCGGCCGCTGGTGGTGTCTATTATCTGGCCTTCGCCGCT  
 TCTGGCTGCCTTACCACCTGGTGAACCTGGTGAAGCCGGCCGCACAGTGGCCGGCTGGGACAAGAACAG  
 CCCTGCGGGCAGCGTCTGAGGCTGGCCGCTATGTGCTCATCGCGTGGCCTTCTGAGCAGCAGCGTG  
 AACCCGGTGTGTACGCGTGGCGGGCGGCCGCTGCTGCGTTCGGCGGGCGTGGGCTTCTGGTCAAGC  
 TACTGGAGGGCACTGGCTCGGAGGTGTCCAGCACCCGCCGGGGCACTCTGGTCCAGACCCCGAAGGA  
 CACACCTGCCTGTCTGAGCCTGGCCCCACCGACAGCTTCATGACTTCTCCACCATTCTGAGTCTTCG  
 AAG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)



<b>ORF Size:</b>	1053 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_008519.2</a> , <a href="#">NP_032545.1</a>
<b>RefSeq Size:</b>	1392 bp
<b>RefSeq ORF:</b>	1056 bp
<b>Locus ID:</b>	16995
<b>UniProt ID:</b>	<a href="#">O88855</a>
<b>Cytogenetics:</b>	14 28.19 cM
<b>Gene Summary:</b>	Receptor for leukotriene B4, a potent chemoattractant involved in inflammation and immune response.[UniProtKB/Swiss-Prot Function]