

Product datasheet for **MG205467**

Aamp (BC031795) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Aamp (BC031795) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Aamp
Synonyms:	MGC30660
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205467 representing BC031795 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGGGCCCGGATGATAGTGAGGTCACCTTTGCATTGCACTCAGCGTCTGTGTTCTGTGTGAGCCTGG
ACCCAAAACCAACACCTTGGCAGTGACAGGGGCGAAGACGACAAGGCCTTCGTGTGGAGGCTCGGTGA
TGGAGAGCTGCTTTGAGTGTGCAGGCCACAAAGACTCTGTGACATGCGCTGGTTTCAGTCACGATTCC
ACCCTAGTGGCCACAGGGGACATGAGTGGTCTTTTGAAGTATGGCAGGTAGACACCAAGGAAGAAGTCT
GGTCTTTGAAGCAGGAGACCTGGAGTGGATGGAGTGGCACCCCTCGGGCTCCTGTCTGTGGCGGGCAC
AGCCGATGGCAACACTTGGATGTGGAAGGTCCCGAATGGTACTGTAAGACCTTTAGGGCCCAACTGC
CCTGCAACCTGTGGCCCGCTCCTCCCTGATGGGAAGAGAGCTGTGGTAGGGTATGAAGATGGTACCATCA
GGATCTGGGACCTCAAGCAAGGAAACCCATCCATGTAATAAAGGGACTGAGGGTCATCAGGGCCCTCT
GACCTGTGTAGCCACCAACCAGGATGGCAGTCTGATCCTGACTGGCTCTGTAGACTGCCAGGCCAAGCTG
GTCAGTGCCACCACTGGCAAGGTGGTGGTGTGTTTCAGACCTGAGACAGTCGCCTCCAGCCAGCCTGG
GAGAGGGGGAGGAGAGTGAAGTCTAACTCCGTGGAGTCTTGGGCTTCTGCAGTGTATGCCTCTGGCCG
TGTTGGCTATCTGGACGGAACCTTGGCCATCTATGACCTGTCTACGCAGACACTCCGGCACCAATGTCAG
CACCAGTCGGGCATTGTTTCAGTCTGTGGGAGGCTGGCACTGCTGTGGTGTACACCTGTAGCCTGGATG
GTGTTGTACGTCTCTGGGATGCCCGACTGGCCGCTGCTTACTGACTACCGGGGCCACACTGCTGAGAT
CCTGGACTTTGCCTCAGCAAAGATGCTTCCCTGGTGGTACCACATCAGGAGACCACAAAGCAAAGTA
TTTTGTGTCCAGAGACTGACCGT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG205467 representing BC031795
 Red=Cloning site Green=Tags(s)

```
MEGPDDSEVTFALHSASVFCVSLDPKTNTLAVTGGEDDKAFVWRLGDGELLFECAGHKDSVTCAGFSHDS
TLVATGDMSSGLLKVWQVDTKEEVWSFEAGDLEWMEWHPRAPVLLAGTADGNTWMWKVPNGDCKTFQGPNC
PATCGRVLDPGKRAVVGIEDGTIRIWDLKQGNPIHVLKGTGEHQGPLTCVATNQDGLILTGSVDCQAKL
VSATTGKVVGVFRPETVASQPSLGESESESNVESLGFCSVMPLAAVGYLDGTLAIYDLSTQTLRHQCQ
HQSGIVQLLWEAGTAVVYTCSLDGVVRLWDARTGRLLTDYRGHTAEILDFALSKDASLVVTTSGDHKAKV
FCVQRPDR
```


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: BC031795

ORF Size:	1076 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC031795 , AAH31795
RefSeq Size:	1535 bp
RefSeq ORF:	1076 bp
Locus ID:	227290
Cytogenetics:	1 38.53 cM