

Product datasheet for **MG205537**

Eif3f (NM_025344) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eif3f (NM_025344) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Eif3f
Synonyms: 0610037M02Rik; AA409853; Eif3s5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG205537 representing NM_025344
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTTCTCCGGCCGTACCGGCTAATGTCCCTCCTGCCACTGCAGCCGAGCCCCGGCGCCGGTTCGTCA
 CCGCAGCCCCGGCTTCAGCCCCGACCCCATCCACGCCAGCTCCGACACCGGCTGCGACTCCCCTGCGTC
 CCCGGCGCCCGTCTCGTCTGATCCTGCTGTAGCTGCGCCTGCAGCCCCGGGCCAGACCCAGCCTCCGCG
 CCAGCCCCAGCGCAGACCGCGCCCTTCGAGCCGGGCCCGCCCTCCCGGGCCCTTCCCGGGCGGCC
 GCGTGGTCAGGCTACACCCCGTCATTTTGGCCTCGATCGTGGACAGCTACGAACCGCGAACGAGGGAGC
 TGCCCCGAGTTATTGGAACCTGTTGGAACTGTTGACAAGCACTCGGTAGAAAGTACCAACTGCTTTTCG
 GTGCCACACAATGAGTCAGAAGATGAAGTGGCTGTTGACATGGAATTTGCTAAGAATATGTATGAATTAC
 AAAAAAAGTCTCCCAAATGAGCTCATCTAGGCTGGTATGCCACAGGCCATGACATCAGAGAACACTC
 AGTGCTGATCCATGAGTACTACAGCAGGGAGGCCCGAACCCATTACCTCACGGTGGACACAGGTCTC
 CAGCATGGGCGCATGAGCATCAAGGCCTATGTCAGCACTTTAATGGGTGTCCCTGGGAGGACCATGGGAG
 TGATGTTACACCTCTCACAGTGAAGTACGCGTATTATGACACTGAACGATTGGAGTTGACCTCATCAT
 GAAGACGTGTTTTAGCCCCAACCGGGTGATTGGACTCTCAAGTGACTTACAACAAGTGGGAGGGCCCTCA
 GCTCGCATCCAGGATGCTAAGCACTGTATTACAGTATGCTGAGGATGTGCTGTGGAAAGTGTCTG
 CTGACAACACGGTGGGCGCTTCTTGATGAGCCTTGCAACCAAGTACCAAGATAGTTCCTGATGACTT
 TGAGACCATGCTCAACAGCAACATCAATGACCTGCTGATGGTGACCTACCTGGCCAATCTCACCCAGTCA
 CAGATTGCCCTCAACGAGAACTTGTAACCTG

AGCGGACCGACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG205537 representing NM_025344
 Red=Cloning site Green=Tags(s)

MASPAVPANVPPATAAAAAPAVVTAAPASAPTSTPAPTAPATPAASPAPVSSDPAVAAPAAPGQTPASA
 PAPAQTPAPSQPGPALPGFPFGRVRLHPVILASIVDSYERRNEGAARVIGTLLGTVDKHSVEVTNCF
 VPHNESEDEVAVDFAKNMYELHKKVSPNELILGWYATGHDITEHSLIHEYREAPNPIHLTVDTGL
 QHGRMSIKAYVSTLMGVPGRTMGVMFTPLTVKYAYDTERIGVDLIMKTCFSPNRVIGLSSDLQQVGGAS
 ARIQDALSTVLQYAEDVLSGKVSADNTVGRFLMSLVNQVPKIVPDDFETMLNSNINDLLMVTYLANLTQS
 QIALNEKLVNL

SGPTRRRLE - GFP Tag - V

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_025344

ORF Size: 1083 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025344.2](#), [NP_079620.2](#)

RefSeq Size: 1252 bp

RefSeq ORF: 1086 bp


Locus ID: 66085

UniProt ID: [Q9DCH4](#)

Cytogenetics: 7 E3

Gene Summary: Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis. The eIF-3 complex associates with the 40S ribosome and facilitates the recruitment of eIF-1, eIF-1A, eIF-2:GTP:methionyl-tRNAⁱ and eIF-5 to form the 43S pre-initiation complex (43S PIC). The eIF-3 complex stimulates mRNA recruitment to the 43S PIC and scanning of the mRNA for AUG recognition. The eIF-3 complex is also required for disassembly and recycling of post-termination ribosomal complexes and subsequently prevents premature joining of the 40S and 60S ribosomal subunits prior to initiation. The eIF-3 complex specifically targets and initiates translation of a subset of mRNAs involved in cell proliferation, including cell cycling, differentiation and apoptosis, and uses different modes of RNA stem-loop binding to exert either translational activation or repression.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG205537