

## Product datasheet for **MG205549**

### Timd4 (BC032879) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Timd4 (BC032879) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Timd4
Synonyms:	Tim4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205549 representing BC032879 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAATTTGATGATTCAAGGCCATCGTTTATTAATAACTAATTACCTCGTGCCGAATTCGGCACGAGGGGC  
TTCTCATCCTCTGGCTGGTGACGGAGCTCTGGTGGCTTTATCTGACACCAGCTGCCTCAGAGGATAACAAT  
AATAGGGTTTTGGGCCAGCCGGTGACTTTGCCTTGTCATTACCTCTCGTGGTCCAGAGCCGCAACAGT  
ATGTGCTGGGGCAAAGGTTTCATGTCCCAATCCAAGTCAATGCAGAGCTTCTCCGTACAGATGGAACAA  
GAATCATCTCCAGGAAGTCAACAAAATATACACTTTTGGGAAGGTCCAGTTTGGTGAAGTGTCTTGAC  
CATCTCAAACACCAATCGAGGTGACAGTGGGGTACTGCTGCCGTATAGAGGTGCCTGGCTGGTTCAAT  
GATGTCAAGAAGAATGTGCGCTTGGAGCTGAGGAGAGCCACAACAACCAAAAAACCAACAACACCCCC  
GGCCAACCACCACCCCTTATGTGACCACCACCACCCAGAGCTGCTTCCAACAACAGTCATGACCACATC  
TGTTCTCCAACCACCACACCACCCAGACACTAGCCACCAGTGCCTTCAAGTACAGCAGTGACCACGTGC  
CCCTCAACAACACCTGGCTCCTTCTACAAGAAACCACAAAAGGGTCCGCCTTCACTACAGAATCAGAAA  
CTCTGCCTGCATCCAATCACTCTCAAAGAAGCATGATGACCATATCTACAGACATAGCCGTACTCAGGCC  
CACAGGCTCTAACCTGGGATTCTCCATCCACTTACAGCTGACGACACAGAAAAACAACATTAACAACA  
AGTGAGCTTTGCAGAAGACAATAAATCACATCAGATCAACAGCAGACAGACCATCTTGATCATTGCCT  
GCTGTGGGATTTGTGCTAATGGTGTTATTGTTTCTGGCGTTTCTCCTTCGAGGGAAAGTCACAGGAGC  
CAACTGTTTGCAGAGACACAAGAGGCCAGACAACACTGAAGATAGTGACAGCGTCTCAATGACATGTCA  
CACGGGAGGGATGATGAAGACGGGATCTCACTCTC

**ACGGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG205549 representing BC032879  
 Red=Cloning site Green=Tags(s)

MNLMIQGHRLKLLITSCRIRHEGLLILWLVTLEWWLYLTPAASEDTIIGFLGQPVTLPCHYLSWSQSRNS  
 MCWKGKSPNSKNAELLRTDGTIRIISRKSTKYTLGKQVQFGEVSLTISNTNRGDSGVYCCRIEVPWFN  
 DVKKNVRLLELRATTTTKPTTTTRPTTTPYVTTTTPELLPTTVMTTSVLPTTTPPQLATTAFASTAVTTC  
 PSTTPGFSQETTKGSFAFTTESETLPASNHSQRSMMTISTDIAVLRPTGNSNPGILPSTSQLTTQKTLTLLT  
 SESLQKTTKSHQINSRQITILIIACCVGFVLMVLLFLAFLLRKGKVTGANCLQRHKRPDNTEDSDSVLNDMS  
 HGRDDEDEGIFTL

TRTRPLE - GFP Tag - V


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

BC032879

<b>ORF Size:</b>	1088 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">BC032879</a>
<b>RefSeq Size:</b>	2331 bp
<b>RefSeq ORF:</b>	1088 bp
<b>Locus ID:</b>	276891
<b>Cytogenetics:</b>	11 B1.1
<b>Gene Summary:</b>	Phosphatidylserine receptor that enhances the engulfment of apoptotic cells. Involved in regulating T-cell proliferation and lymphotoxin signaling. Ligand for HAVCR1/TIMD1. [UniProtKB/Swiss-Prot Function]