

Product datasheet for **MG205657**

Rbm39 (BC030493) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rbm39 (BC030493) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rbm39
Synonyms:	1500012C14Rik; 2310040E03Rik; B330012G18Rik; C79248; caper; R75070; Rnpc2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205657 representing BC030493 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAACTGGCAGCAAGAATTCGACCGAGAGATTTGGAAGAGTTTTCTCTACAGTTGGAAGGTTTCGAG
ATGTGAGAATGATTTCTGATAGAAATTCAGACGTTCCAAAGGAATTGCATATGTGGAATTTGTTGATGT
CAGTTCAGTGCCTCTAGCAATTGGATTAAGTGGCCAACGAGTTTTAGGAGTGCCAATCATAGTTCAGGCA
TCACAGGCAGAAAAACAGAGCTGCTGCAATGGCAAACAATTTACAAAAGGGAAGTCTGGACCTATGA
GGCTCTATGTAGGCTCATTACACTTTAACATAACTGAAGACATGCTTCGGGGGATCTTTGACCTTTTGG
AAGGATTGAAAGTATTCAGCTCATGATGGATAGTAAACTGGCCGATCTAAGGGATATGGATTTATTACA
TTTTCTGACTCAGAATGTGCCAAAAAGGCTTTGGAACAACCTGAATGGATTTGAGTTAGCTGGGAGACCAA
TGAAAGTTGGTCACGTTACTGAACGAACTGATGCTTCGAGTGCTAGTTTCATTTTTGGACAGTGACGAACT
AGAAAGGACTGGAATTGACTTGGGAACAACCTGGACGTCTTCAGTTAATGGCAAGACTTGGCGAGGGTACC
GGTTTGCAGATTCACCTGCCGCACAACAAGCATTACAAATGAGTGGCTCTCTGGCATTGGTGTGTGG
CAGATTTGCAAACACGACTTTCCAACAGACCGAAGCCTCAGCTTAGCTGCAGTGCCTGTGTTCAACC
TCTTGAACACAGTGTTCCTCAACTTTCTAACATGTTAATCCTCAAACAGAAGAAGAAGTTGGATGGGAT
ACAGAGATTAAGATGATGTCATTGAAGAATGTAATAAACATGGAGGAGTCATTCATATTTATGTTGATA
AAAATTCAGCTCAGGGCAATGTGTATGTGAAGTGCCCATCTATTGCTGCGGCTATTGCTGCTGCTCAATGC
ATTGCATGGCAGATGGTTTGTGGTAAATGATAACAGCAGCGTATGTACCTCTTCCAACCTACCACAAC
CTCTTCTGATTCTATGACAGCAACAACACTACTGGTTCCAAGTAGACGA

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG205657 representing BC030493
 Red=Cloning site Green=Tags(s)


MQLAARIRPRDLEEFFSTVGKVRDVRMISDRNSRRSCKGIAYVEFVDVSSVPLAIGLTGQVRLGVPIIVQA
 SQAENRRAAMANNLQKGSAGPMRLYVGSLSHFNITEDMLRGI FEPFGRIESIQLMMDSETGRSKGYGFI T
 FSDSECAKKALEQLNGFELAGRPMKVGHVHTERTDASSASSFLDSELERGTGIDLTTGRLQLMARLAEGT
 GLQIPPAQQALQMSGSLAFGAVADLQTRLSSQTEASALAAAASVQPLATQCFQLSNMFNPQTEEEVGVWD
 TEIKDDVIEECNKHGGVITHIYVDKNSAQNVVYVKCP SIAAAIAAVNALHGRWFAGKMITAAYVPLPTYHN
 LFPDSMTATQLLVPSRR

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

BC030493

ORF Size:	1103 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC030493 , AAH30493
RefSeq Size:	2413 bp
RefSeq ORF:	1103 bp
Locus ID:	170791
Cytogenetics:	2 H1
Gene Summary:	Transcriptional coactivator for steroid nuclear receptors ESR1/ER-alpha and ESR2/ER-beta, and JUN/AP-1. May be involved in pre-mRNA splicing process.[UniProtKB/Swiss-Prot Function]