

Product datasheet for **MG205724**

Mtch1 (BC023135) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mtch1 (BC023135) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mtch1
Synonyms:	C77849, PSAP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205724 representing BC023135 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGGAGCCGGAGCTGGTGCAGGAGCTCGGGCGGCGCGCCGGCGGAGTCGAGGCCCGCGCTCGGG
ACCCACCAACCGCGCACCGCGCACCCCTCGCCATCCTCGGCCCGCGGCTCAGCCGTCGGCGCGCAGGAT
GGACGGCGGCCCGGGCGCCCGGGCTCCGGGACAACGCCCGACCAACCGAGGCGCTGTTCTGTGGCGCTG
GGCGGGCGGTGACGGCTCTCAGTCACCGCTGCTCTACGTGAAGCTGCTGATCCAGGTGGGTGATGAGC
CGATGCCCCCACCTTGGGACCAATGTGCTGGGAGGAAGGTCTCTACCTGCCGAGCTTCTTCACCTA
TGCCAAGTACATTGTGCAGGTGGATGGGAAGATAGGGCTCTTCCGGGGCTGAGCCCCCGCTTATGTCC
AACGCCTTGTCCTACTGTGACCCGCGGCAGCATGAAGAAGTTTTCCCTCCAGATGAGATGGAGCAGGTTT
CCAACAAGGACGACATGAAGACCTCACTCAAGAAAGTTGTGAAGGAGACATCGTATGAGATGATGATGCA
GTGTGTATCGGAATGCTGGCCCATCCCTACACGTGATCTCGATGCGATGCATGGTGCAGTTTGTGGGA
CGGAGGCCAAGTACAGTGGTGTGCTGAGTTCTATTGGGAAGATCTTCAAGGAAGAGGGGCTGCTGGGAT
TCTTCGTTGGCTTAATCCCTCACCTCCTGGCGATGTGGTTTTCTGTGGGGCTGTAACCTGCTGGCCCA
CTTCATCAATGCCTACTTGGTGGACGACAGCGTGAGTGACACCCAGGGGGCTGGGAAACGACCAGAAT
CCAGGTTCCAGTTTAGCCAGGCCCTGGCCATCCGGAGCTACACCAAGTTTGTGATGGGATTGCACTGA
GCATGCTGACCTACCCCTTCTGCTCGTTGGAGATCTCATGGCAGTGAACAACTGTGGGCTGCGGGCTGG
ACTCCCTCCGTATTCCCCTGTGTTCAAGTCCTGGATCCACTGCTGGAAGTACCTGAGTGTGACGGGCCAG
CTCTCCGCGGCTCCAGCCTGCTTTCCGCCGGGTGTCATCGGGGTGATGCTTTGCCCTGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG205724 representing BC023135
 Red=Cloning site Green=Tags(s)

MAGAGAGAGARGGAPAGVEARARDPPPAHRAHPRHPRPAAQPSARRMDGGPGAPGSGDNAPTTEALFVAL
 GAGVTALSHPLL YVKLLIQVGHEPMPPTLGTNLGRKVL YLP SFFTYAKYIVQVDGKIGLFRGLSPRLMS
 NALSTVTRGSMKKVFPDEMEQVSNKDDMKTSKKVVKETSYEMMMQCVSRMLAHPLHVISMRCMVQFVG
 REAKYSGVLSSIGKIFKEEGLLGFFVGLIPHLLGDVFLWGCNLLAHFINAYLVDDSVSDTPGGLGNDQN
 PGSQFSQALAIRSYTKFVMGIAVSMLTYPFLLVGDLMAVNNCGLRAGLPPYSPVFKSWIHCWKYLSVQGG
 LFRGSSLLFRRVSSGSCFALE

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC023135

ORF Size: 1113 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC023135.1](#)

RefSeq Size: 1776 bp

RefSeq ORF: 1115 bp

Locus ID: 56462

Cytogenetics: 17 A3.3

Gene Summary: Potential mitochondrial transporter. May play a role in apoptosis (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG205724