

Product datasheet for **MG205986**

Col9a3 (BC030945) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Col9a3 (BC030945) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Col9a3
Synonyms:	AV006866
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG205986 representing BC030945 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAACTACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGCCAGGCAAGGATGGCAAAGATGGTGTGCCGGGACTTGATGGTGAGAAGGGAGAGGCTGGTCGCAATG
GTGGCCAAAGGAGAGAAAGGCCCAATGGGCTGCCGGGGCTCCCTGGACGAGCAGGGTCCAAGGCGAGAA
GGGAGAACCGGGTAGAACTGGCGAGCTGGGTGAGGCTGGCCCTCTGGAGAGCCAGGTATCCCTGGAGAT
GTTGGTGTCCAGGGGAGCGTGGTGAGGCTGGTCACAGGGGCTCCGTGGGAGCTCTTGCCCAACAAGGTC
CTCCTGGTGCTCCTGGCATCCGTGGCTTTCAGGGACAAAAAGGCAGCACAGGAGACCCTGGCCTTCCGGG
CCCCCAAGGCCTACGGGGAGATGTGGGTGACCGGGGTCCAGGAGGTGCCACAGGCCCTAAGGGAGACCAG
GGCATTGCAGTTCCGATGGTCTTCCAGGGGACAAAGGCAGCTGGGTCTAATGGCCCCGTTGGACAAA
AAGGAGAGTCTGGCAGCCGAGGGGAAGTGGGCCCAAGGCATCCAGGGCCCCAATGGCACCAGTGGAGT
CCAGGGGTGTACCTGGTCCCCAGGTCCACTGGGCCTCCAAGGTGTGCAGGGTGTCCAGGCATCACCGGG
AAGCCTGGAGTTCGGGCAAGGAAGCCAGTGAACAGCGCATCAGGGAGCTATCGGGGGTATGATCAGTG
AGCAAATTGCACAGTTGGCTGCACACCTGAGGAAGCCCTAGCACCAGGCTCCATCGGCAGGCCTGGTCC
AGCTGGCCCCCAGGCCCTCAGGCCCTCAGGCTCTATTGGCCACCCTGGCGCTCGGGTCCCCCTGGA
TACCGTGGTCCCACTGGGGAAGTGGAGATCCTGGACCAAGAGGGAGCCAGGGTGACAGAGGGGACAAGG
GAGCAACAGGTGCAGGGCTGGATGGGCCTGCTGGGGACCAGGGCTACCAAGGGCCTCAAGGTGTACCTGG
CATCAGCAAAGATGGTCGAGATGGTGCTCATGGTGAGCCTGGCCTTCCCGGTGATCCTGGCCTTCTGGA
GCTGCTGGTGCTCAGGGTACCCAGGGATCTGCGACACTTCAGCCTGCCAAGGAGCTGTGTAGGGGGAG
GTGGGGAAAAGTCGGGTCTCGGAGCTCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >MG205986 representing BC030945
 Red=Cloning site Green=Tags(s)


MPGKDGKDGVPGLDGEKGEAGRNGGQGEKGNLPLPLPGRAGSKGEKGEPRGTGELGEAGPSGEPGIPGD
 VGVPPERGEAGHRGSVGLGPQPPGAPGIRGFQKQKSTGDPGLPGPQGLRGDVGDRPGGATGPKGDQ
 GIAGSDGLPGDKGELGPNPVGQKGESESRGELGPKGIQPNGTSGVQGVPPGPGPLGLQGVQGVPGITG
 KPGVPGKEASEQRIRELCGMISEQIAQLAAHLRKPLAPGSIGRPGPAGPPGPPGPPGSIHGPARGPPG
 YRGPTGELGDPGPRGSQGDGRDKGATGAGLDGPAGDQGYQGPQGVPGISKDGRDGAHGEPGLPGDPLPG
 AAGAQTGTPGICDTSACQGAVLGGGGEKSGPRSS

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: BC030945

ORF Size: 1149 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: BC030945, AAH30945


RefSeq Size: 1605 bp

RefSeq ORF: 1151 bp

Locus ID: 12841

Cytogenetics: 2 103.14 cM

Product images:



Circular map for MG205986