

Product datasheet for **MG206083**

Nr2f6 (BC008138) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nr2f6 (BC008138) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Nr2f6
Synonyms:	EAR2, COUP-TF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206083 representing BC008138 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCATGGTGACCGGTGGCTGGGGCGACCCCGAGGCGACACGAACGGCGTGGACAAGGCTGGTGGGA
GCTACCCACGCGCGACCGAGGACGATTCGGCGTCACCTCCCGGGCGACACGCGCGGAGCCGGGCGA
CGAGGAGCGTCCGGGTTGCAGGTGGACTGCGTGGTGTGCGGGGACAAGTCCAGTGGAAAGCATTACGGC
GTGTTACCTGCGAGGGCTGCAAGAGTTTCTTCAAGCGCAGCATCCGCCAATCTCAGTACACCTGCC
GGTCCAACCGTGACTGTCAGATTGATCAGCACCACCGAACCAGTGTGACTGTGCGGCTCAAGAAGTG
CTCCGGGTGGGCATGCGCAAGGAGGCCGTGACGCGAGGCCGCATCCCGCATGCGCTCCCGGTCCAGCG
GCCTGCAGTCCCCGGGCGCGACGGGCGTCAACCTTTCACGGGGCCGCCAGTGTCCGAGCTGATTGCGC
AGCTGCTGCGTGTGAGCCCTACCCCGCGGCCGACGCTTTGGTGGCGGGCGCTGACTGGGCATCGA
CAACGTGTGCGAGTTGGCGGCACGCCTGCTGTTTCAAGCAGGTCGAGTGGGCCCGCCACGCGCCCTTCTT
CCCGAGCTGCCGGCCCGCACCAGGTGGCGTGTGCGGCTCAGTGGAGTGAAGTCTTCGCTGCTGAACG
CGGCGCAGGCGCGCTGCCGCTGCATACGGCACCCTGCTGGCCGCCCGGGGTTGCATGCCCGCCCAT
GGCAGCCGAGCGGGCCGTGGCTTCATGGACCAGTGCCTGCTTCCAGGAGCAGTGGACAAGCTGGGC
CGCTTCAGGTGGATGCTGCGGAGTACGGCTGCCTCAAGGCCATCGCGCTTTCACGCTGATGCCTGTG
GCCTTTCTGACCCAGCCCATGTGGAGAGCCTGCAGGAGAAGGCACAGGTGGCCCTACCGAGTATGTGCG
TGCCAGTACCCATCGCAGCCCCAGCGCTTTGGCGTCTGCTGCTGCGGCTGCCAGCCCTGCGTGTGTG
CCCGCATCCCTCATCTCCAGCTCTTTCATGCGCCTGGTGGCAAGACACCCATCGAGACCCTCATCC
GGGACATGCTTCTGTCAGGGAGCACCTTAACTGGCCCTATGGCTCGGGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG206083 representing BC008138
Red=Cloning site Green=Tags(s)

MAMVTGGWGDPPGGDTNGVDKAGGSYPRATEDDSASPPGATSDAEPGDEERPGLQVDCVVCVCGDKSSGKHYG
VFTCEGCKSFFKRSIRRNL SYTCRSNRDCQIDQHHRNQCQYCRLKKCFRVGMKEAVQRGRIPHALPGPA
ACSPPGATGVEPFTGPPVSELIAQLLRAEPYPAAGRFGGGGAVLGDNVCELAARLLFSTVEWARHAPFF
PELPAADQVALLRLSWSELFVLNAAQAALPLHTAPLLAAAGLHAAPMAAERAVAFMDQVRAFQEQVDKLG
RLQYDAAEYGLKAIALFTP DACGLSDPAHVESLQEK AQVALTEYVRAQYPSQPQRFGRLLLRLPALRAV
PASLISQLFFMRLVGKTP IETLIRDMLLSGSTFNWPYGS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC008138

ORF Size: 1170 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC008138](#), [AAH08138](#)

RefSeq Size: 1959 bp


RefSeq ORF: 1172 bp

Locus ID: 13864

Cytogenetics: 8 34.43 cM

Gene Summary: Transcription factor predominantly involved in transcriptional repression. Binds to promoter/enhancer response elements that contain the imperfect 5'-AGGTCA-3' direct or inverted repeats with various spacings which are also recognized by other nuclear hormone receptors. Involved in modulation of hormonal responses. Represses transcriptional activity of the lutropin-choriogonadotropic hormone receptor/LHCGR gene, the renin/REN gene and the oxytocin-neurophysin/OXT gene. Represses the triiodothyronine-dependent and -independent transcriptional activity of the thyroid hormone receptor gene in a cell type-specific manner. The corepressing function towards thyroid hormone receptor beta/THRB involves at least in part the inhibition of THRB binding to triiodothyronine response elements (TREs) by NR2F6. Inhibits NFATC transcription factor DNA binding and subsequently its transcriptional activity. Acts as transcriptional repressor of IL-17 expression in Th-17 differentiated CD4(+) T cells and may be involved in induction and/or maintenance of peripheral immunological tolerance and autoimmunity. Involved in development of forebrain circadian clock; is required early in the development of the locus coeruleus (LC). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG206083