

## Product datasheet for **MG206087**

### **Ncf1 (NM\_010876) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ncf1 (NM\_010876) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ncf1  
**Synonyms:** Ncf-1; NCF-47K; NOXO2; p47-phox; p47; p47phox  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG206087 representing NM\_010876  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGGGGACACCTTCATTCGCCATATTGCCCTGCTGGGCTTCGAGAAGCGCTTCATCCCCAGCCAGCACT  
 ATGTGTACATGTTCTGGTTAAGTGGCAGGACCTGTCGGAGAAGGTGGTCTACAGAAAATTCACCGAGAT  
 CTACGAGTTCCATAAAAATGCTGAAGGAGATGTTCCCCATTGAGGCCGGCGAGATCCACACAGAGAACAGA  
 GTCATCCACACCTCCCGCACCCAGGTGGTTTGGATGGGCAACGAGCCGCTGAGAGTCGCCAGGGCACTC  
 TCACTGAATACTCAACGGCCTCATGGGACTGCCCGTGAAGATCTCCCGTGCCACACCTGCTGGACTT  
 CTCAAAGTGCGGCCTGATGACCTGAAACTGCCCACTGACAGCCAGGCGAAGAAGCTGAGACATACCTG  
 GTGCCAAAGATGGCAAGAATAACGTAGCTGACATCACAGGCCCATCATCCTTCAGACCTATCGGGCCA  
 TTGCTGACTACGAGAAGAGTTCGGGAACAGAGATGACCGTGGCAACCGGAGACGTGGTGGACGTTGTGGA  
 GAAGAGCGAGAGCGGCTGGTGGTTTTGCCAGATGAAGACAAAGCGAGGTTGGGTCCCTGCATCCTATCTG  
 GAGCCCTTGACAGTCCCGACGAGGCGGAGGATCCGGATCCCAACTACGCAGGTGAACCGTATGTAACCA  
 TCAAAGCGTACGCTGCTGTTGAAGAGGACGAGATGTCCCTGTCTGAGGGTGAAGCCATTGAGGTCATTCA  
 TAAGCTCCTGGATGGCTGGTGGTGGTCAGGAAAGGGGATATCACCGCTATTTCCCATCCATGTATCTG  
 CAGAAGGCTGGGAGGAGATAACCCAGGCCAGCGACAGATCAGAGGCCGCGGGGCCACCTCGAAGGT  
 CGACCATCCGCAACGCACAGAGCATCCACCAGCGTTCTCGGAAGCGCCTTAGCCAGGACACCTATCGCCG  
 CAACAGCGTCCGATTCTGCAGCAGCGCAGACGCCCGGGCGACCCGGGCGCTGAGCACGGATGGCACA  
 AAGGACAATCCATCGACTCCGCGCTCAAACCACAGCCCGCGGTGCCTCCGCGACCCAGCTCGGACCTCA  
 TCCTGCACCGCTGCACAGAGAGCACAAACGGAAGCTGACGTCCGCTGTG

**ACCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG206087 representing NM\_010876  
Red=Cloning site Green=Tags(s)

MGDTFIRHIALLGFEKRFIPSQHYVYMFVLVKWQDLSEKVYRKFTEIYEFHKMLKEMFPIEAGEIHTENR  
 VIPHLPPAPRWFQRAAESRQGLTEYFNGLMGLPVKISRCPHLLDFKVRPDDLKLPDTSQAKKPETYL  
 VPKDGNVADITGPIILQTYRAIADYEKSSGTEMTVATGDVVDVVEKSESGWWFCQMKTKRGWVPASYL  
 EPLDSPDEAEDPDPNYAGEPYVTIKAYA AVEEDEMSEGEAIEVIHKLLDGGWVVRKGDITGYFPSMYL  
 QKAGEEITQAQRQIRGRGAPRRRSTIRNAQSIHQRSRKLRSQDTYRRNSVRFLLQRRRRPGRGPLSTDGT  
 KDNPSTPRVKPQPAVPPRPSSDLILHRCTESTKRKL TSAV

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_010876


**ORF Size:** 1170 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_010876.4</a></u>
<b>RefSeq Size:</b>	2721 bp
<b>RefSeq ORF:</b>	1173 bp
<b>Locus ID:</b>	17969
<b>UniProt ID:</b>	<u><a href="#">Q09014</a></u>
<b>Cytogenetics:</b>	5 74.47 cM
<b>Gene Summary:</b>	NCF2, NCF1, and a membrane bound cytochrome b558 are required for activation of the latent NADPH oxidase (necessary for superoxide production).[UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for MG206087