

Product datasheet for **MG206346**

Pcyt2 (NM_024229) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pcyt2 (NM_024229) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pcyt2
Synonyms:	1110033E03Rik; ET
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206346 representing NM_024229 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGATCCGGAACGGGCACGGGGCTGCCAGCGCTGCTGGGCTGAAGGGTCCGGGAGATCAGCGCATCGTGC
GGGTGTGGTGGCATGGCTGCTATGACATGGTGCATTATGGCCACTCCAACCAGCTACGGCAGGCACGGGC
CATGGGGACTACCTCATCGTGGGTGTGCATACGGATGAGGAGATTGCCAAGCATAAGGGGCCCCCGGTG
TTTACCCAGGAGGAGAGGTACAAGATGGTACAGGCCATCAAGTGGGTGGATGAGGTGGTCCCGCTGCTC
CCTACGTCACCACCTGGAGACTGGACAAGCACAACCTGTGACTTCTGTGTTACAGGCAATGACATCAC
GCTGACAGTAGACGGCCGAGACCTATGAGGAAGTGAAGCAGGCTGGGAGGTACAGAGAGTGCAAACGC
ACCCAGGGTGTGTCCACCACAGACCTCGTGGTGCATGCTGCTGGTGACCAAGGCCATCACAGCAGCC
AGGAGATGTCTCTGAGTACCGGAATATGCTGACAGTTTTGAAAGCCCTCACCCGACACCTGCCGG
GGACACACTTTCCTCAGAAGTCTCTCCAGTGCCTGGGGGCAGAGCCCTGGACAGGGGTGTCCAG
TTTCTACAGACATCCAGAAGATCATCCAGTTTGTCTTGGGAAGGAGCCCGCCGGGAGACGGTCA
TCTATGTGGTGGTGCCTTTGACCTGTTCCACATCGGGCAGCTGGACTTCTACAGGAGGTGCACAAGCT
AGCCAAGAGGCCCTACGTCATCGCCGGCCTACACTTTGACCAGGAAGTAAACCGGTACAAGGCAAGAAC
TACCCATCATGAACCTGCACGAGCGGACTCTCAGTGTGCTGGCCTGCCGTATGTTTCAGAAGTGGTGA
TTGGGGCACCATACTCCGTGACAGCGGAGCTCCTGAATCACTTCAAGGTGGACCTGGTGTGTCACGGGAA
GACAGAGATTGTACCCGACAGGGATGGGTCTGACCCCTACCAGGAGCCCAAGAGAAGAGGATCTTCTAT
CAGATTGACAGTGGCAGTGACCTCACTACAGACCTGATTGTGCAGAGGATCATCAAGAACAGGCTGGAGT
ATGAAGCACGGAATCAGAAGAAAGAAGCCAAGGAATTAGCCTTTCTGGAGGCCACGAAGCAGCAGGAGGC
GCCGCTGGAGGGGAGATTGAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >MG206346 representing NM_024229
 Red=Cloning site Green=Tags(s)

MIRNGHGAASAAGLKGPGRIVRVWCDGCDMVIHYGHSNQLRQARAMGDYLIVGVHTDEEIAKHKGPPV
 FTQEERYKMVQAIKWVDEVVPAAPYVTTLETLDKHNCDFCVHGNDITLTVDGRDYEYEVKQAGRYRECKR
 TQGVSTTDLVGRMLLVTKAHHSSQEMSSEYREYADSFQKPPHPTPAGDTLSSEVSSQCPGGQSPWTGVSQ
 FLQTSQKIIQFASGKEPQPGETVIYVAGAFDLFHIGHVDFLQEVHKLAKRPYVIAGLHFDQEVNRYKGN
 YPIMNLHERTSLVLAARYVSEVVIYVAGAFDLFHIGHVDFLQEVHKLAKRPYVIAGLHFDQEVNRYKGN
 QIDSGSDLTTDLIVQRIKRNRYEARNQKKEAKELAFLEATKQEQEAPPGGEID


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_024229

ORF Size:	1212 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024229.3
RefSeq Size:	1908 bp
RefSeq ORF:	1215 bp
Locus ID:	68671
UniProt ID:	Q922E4
Cytogenetics:	11 E2
Gene Summary:	Plays an important role in the biosynthesis of the phospholipid phosphatidylethanolamine. Catalyzes the formation of CDP-ethanolamine (By similarity).[UniProtKB/Swiss-Prot Function]