

Product datasheet for **MG206469**

1810007M14Rik (BC014838) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	1810007M14Rik (BC014838) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	1810007M14Rik
Synonyms:	1810007M14Rik; C21orf66; Gcfc; Gcfc1; Pax3/7bp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206469 representing BC014838 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCACCAATCTTGATTCCTTTGGACGTGATCGGGCACTTTATCAAGAACATGCAAAACGTCGGATTG
CAGAGCGGAAGCCAGGAGGACTCGCCGTAGGCAAGCCAGAGAGCAAACCGGTCAGATGGCAGATCACCT
GGAAGGCTGTCCAGTGATGATGAAGAGACGTCTACAGATATTACTAATTTCAATCTGGAGAAAGATCGA
ATTTTGAAGAGTCTAGCAAAGTTTTGAAGATGTCCTTGAGAGTTTCTATTCAATTGACTGCATTAAG
CACAGTTTGAAGCATGGCGTTCAAATACTACATGTCCTATAAGGATGCTTATATTGGTCTTTGTCTGCC
CAAAGTGTCAACCCCTAATCAGACTTCAGCTTCTCACCTGGACTCCTCTTGAGGCAAAATGCCGTGAC
TTTGAAACCATGCTGTGTTTGAATCTCTGCTGTTCTATGGCTGTGAAGACCGAGAGCAGGAGAAGGATG
AGGCTGACGTTGCTCTGCTGCCACCATTTGTTGGAAGGTCATTCTCCCAAGCTGACAGTGATAGCTGA
AACCATGTGGGATCCCTTTCTACAACACAGACTTCAAGGATGGTTGGAATCACAATGAAATTAATAAT
GGCTACCCCTCAGTAGTCAATGCAGACAATAAAACACGCAGGTATACCTTAAAGCTCTTTACTGAGGA
TGAGAAGAACTTTAGATGATGATGATTTATGCCCTATATCCCAAAATGTCTTAGAAAATAAAATTC
TGGCCTTACTTATTTTTCAACGACAGTTTTGGTCTGCTGTTAAGCTATTGGGAAATTTCTTCAGTGG
TATGGCATTTTCTCAACAAAACCTTTCCAGGAGTTATCAATAGATGGCTTGTTAAACAGATACATCCTTA
TGGCTTTTCAGAACTCAGAATACGGAGATGACAGCATCAGAAAAGCCCAAAATGTAATCAACTGTTTCCC
TAAACAGTGTTTGTGAATCTGAAAGGAGAAAGAACTATTTCTCAGTTAGAAAATTTTGGCGATATCTT
GTACATTTAGCAGATACAATTTACAGGAACAGTATTGGGTGCTCTGATGTGGAAGAAAGAAATGCAAGGG
AAAATATAAACAGATAGTAAACTTCTTGCAAGTGTTCCGGCTTTGGATCACGCTATATCTGTTGCAAG
TGACCACAATGTGAAGGAAGTCAAGTCTCTGATTGAAGGAAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)

Protein Sequence: >MG206469 representing BC014838
 Red=Cloning site Green=Tags(s)


MAPNLDSFGRDRALYQEHAKRIAEREARRTRRRQAREQTGMADHLEGLSSDDEETSTDITNFNLEKDR
 ILKESSKVFEDVLESFYSIDCIKAQFEAWRSKYYSYKDAYIGLCLPKLFNPLIRQLLTWTPLEAKCRD
 FETMLWFESLLFYGCEDREQEKDEADVALLPTIVEKVILPKLTVIAETMWDPFSTTQTSRMVGITMKLIN
 GYPSVYNADNKNTQVYLKALLLRMRRTLDDVFMPLYPKNVLENKNSGPYLFFQRQFWSSVKLLGNFLQW
 YGIFSNKTLQELSIDGLLNRYILMAFQNSEYGGDSIRKAQNVINCFPKQWFVNLKGERTISQLENFCRYL
 VHLADTIYRNSIGCSDVEKRNARENKIQIVKLLASVRALDHAISVADSHNVKEVKSLEIGK

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: BC014838

ORF Size: 1233 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC014838.1](#)

RefSeq Size: 1414 bp


RefSeq ORF: 1235 bp

Locus ID: 67367

Cytogenetics: 16 C3.3

Gene Summary: Adapter protein linking the transcription factors PAX3 and PAX7 to the histone methylation machinery and involved in myogenesis. Associates with a histone methyltransferase complex that specifically mediates dimethylation and trimethylation of 'Lys-4' of histone H3. Mediates the recruitment of that complex to the transcription factors PAX3 and PAX7 on chromatin to regulate the expression of genes involved in muscle progenitor cells proliferation including ID3 and CDC20.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG206469