

Protein Sequence: >MG206474 representing NM_019794
 Red=Cloning site Green=Tags(s)

MANVADTKLYDILGVPPGASENELKKAYRKLAKKEYHPDKPNAGDKFKEISFAYEVL SNPEKREL YDRYG
 EQGLREGSGGGGGMDDIFSHIFGGGLFGFMGNQSRSRNGRRRGGEDMMHPLKVSLEDLYNGKTTKLQLSKN
 VLCSACSGQGKSGAVQKCSACRGRGVRIMIRQLAPGMVQQMQSVCSDCNGEGEVINEKDRCKKCEGKKV
 IKEVKILEVHVVDKGMKHGQRITFTGEADQAPGVEPGDIVLLQLQEKEHEVFQRDGNLHMTYKIGLVEALC
 GFQFTFKHL DARQIVVKYPPGKVI EPGCVRVVRGEGMPQYRNPFKGDLYIKFDVQFPENNWINPKLSE
 LEDLLPSRPEVPNVIGETEEVELQFEDSTRSGGGQRREAYNDSSDEESSHHGPGVQCAHQ

TRTRPLE - GFP Tag - V


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_019794

ORF Size:	1236 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019794.1 , NP_062768.1
RefSeq Size:	2838 bp
RefSeq ORF:	1239 bp
Locus ID:	56445
UniProt ID:	Q9QYJ0
Cytogenetics:	8 C3
Gene Summary:	Co-chaperone of Hsc70. Stimulates ATP hydrolysis and the folding of unfolded proteins mediated by HSPA1A/B (in vitro).[UniProtKB/Swiss-Prot Function]