

Product datasheet for **MG206518**

Snap47 (NM_144521) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Snap47 (NM_144521) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Snap47
Synonyms: 1110031B06Rik; SNAP-47
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG206518 representing NM_144521
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGTTCGATATGCGTGTCCACTCCTGGTCTTGTTCCTACTATCTGGACCTTGAGAAGCAATGGGTTT
 CTGGAAAACACTGACCCTGACTCCTCATTCACTGAAGTTCATTGTGGAAAAACAGAAGAGGTCCTTGTGG
 CCTGCCCTCTCCAGCATCATTGAGATCAGGAAAGAGTCATCTCTTTATCTTCGGTGGCATCACAGTC
 CTGGAGAAGGGCCAGACCAAGCACTGGTTCAGCTCACTGCAGCCTAGTCGAAATGTGGTCTTCAACGTCA
 TCGAACATTTCTGGAGGGAGCTGCTACTGTCCCAGCCAGGCACTGCTGCCAACATACCCTCCCACGTGAC
 CAGGGGACAGGAGCTGATAGGACTGATGGCAAATCCCAGAAACGCATGGAAGACACAGCGAAGGATCTC
 CAGCAGCAGAGCGAACAGCTGGACAGCGTCTGAAGGGCCTGGAGAAGATGGAGTCCGACCTGGACGTGG
 CTGACAGATTGCTGACAGAATTGGAGACACCTTCTGGTGGCCTTTTGGCTCTAAGTTTTGGAAGATGCC
 GGCAGAAGAAAATCTCAAGGAAGGTGTGTCTTCCACTTGTGAGCCCTTTGGAAAGGAAGGGGTTGTGATC
 ACGGTCCCTGCTATTATTTCTGAGAGAGCAGAGTCTATTCCAAGCTGGGAAAGCTCACTGTCTGGTCT
 CCGCTCTGGAAATCTATGACTCCTGCTCATTGCTCCTGCACAGATTTGAGAAGGAAGATGTTGATGATAT
 CAAGGTTCACTCTCCTTATGAGGTCAAGTACGCCAGCGCTTATTGGGAAGCCAGACGTGGCTTACCAG
 CTGATATCTGCCAAGATGCCGGAGGTCAATCCCATCTTAGAAGTACAGTTCAGCAGCAAGATTGAGCTCC
 TGGAAGACGCCTTAGTGCTCAGGAACAAGGTCTTTGCTTCTCTGCAGAGAGGCATGCAGCATCCAGACC
 GAAGGGCTGCACTCCACACCGAGAGCTGCCACAGGAGGACAGGAAGGCGAGCAGCTGCAGCTGCAGAAG
 AACTTGCCACTGTTCTCTGAGGGGGAAGCCAGGAGCTGACACAGATCCTGAGCAAAATGAAGGGACTCG
 CACTGGACACTGAGGCAGAGCTGGAACGCCAGGATGCGGCTCTGGATGGCATCACCGTAGCTGTGGACCC
 GGCAACCTGAATGTGGACAAGCAAAACCGGCGTATGAGAAAGCTGATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG206518 representing NM_144521
Red=Cloning site Green=Tags(s)


MSSDMRVHSWCSYYLDLEKQWVSGKLTLPKSLKFLVEKTEEVLVGLPLSSIIIEIRKESLFFIGAITV
 LEKGQTKHWFSSSLQPSRNVVFNVEHFWRLLLSQPGTAAIPSHVTRGQELIGLMANSQKRMEDTAKDL
 QQQSEQLDSVLKGLEKMESDLVDADRLLELETPSWWPFSGKFWKMPAEENLKEGVSSTCEPFKGEVVI
 TVPAIISERAESHKLGKLTVLVSALEIYDSCSLLLHRFEKEDVDDIKVHSPYEVSIQRQFIGKPDVAYQ
 LISAKMPEVPILEVQFSSKIELLEDALVLRNKVFASSAERHAASRPKGCTPHRELPTGGQEGEQLQLQK
 NLPLFSEGEAQELTQILSKMKGLALDTEAELERQDAALDGITVAVDVDRATLNVDKQNRMRKLM

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja3476_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_144521

ORF Size: 1239 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144521.2](#), [NP_653104.1](#)

RefSeq Size: 1837 bp

RefSeq ORF: 1242 bp


Locus ID: 67826

UniProt ID: [Q8R570](#)

Cytogenetics: 11 B1.3

Gene Summary: May play a role in intracellular membrane fusion.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG206518