

Product datasheet for **MG206549**

Ssb (NM_009278) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ssb (NM_009278) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ssb
Synonyms:	SS-B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206549 representing NM_009278 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGAAAAATGGAGATAATGAAAAATGACTGCTTTGGAGGCCAAAAATCTGTCATCAAATTGAGTATT
ATTTTGGAGACTTCAATTTGCCACGAGACAAGTTTTTAAAAGAACAGATCAAATTGGATGAAGGCTGGGT
ACCTTTGGAAACAATGATAAAATCAACAGGCTAAACCGCTGACAACAGACTTTAATGTAATTGTGCAA
GCACTGAGCAAATCTAAGGCAAACTCATGGAAGTCAGTGCAGACAAAATAAAATTAGAAGATCACAA
GCAGACCACTCCCTGAAGTGACGGATGAGTATAAGAATGATGTA AAAAACAGATCTGTTTATATTAAGG
TTTCCCAACTGACGCCACCTTGATGATATAAAAGAATGGCTAGACGATAAAGGCCAAATACTGAATATT
CAAATGAGAAGAACATTACAAAAACATTTAAGGGTCAATATTTGCTGTGTTTGATAGTATTCAGTCTG
CAAAGAAGTTTGTGGAGATCCCTGGCCAGAAGTACAAAGACACTAACCTGCTAATACTCTTTAAGGAAGA
TTACTTTGCAAAAAAAAAATGAAGAAAGAAAGCAGAGCAAAGTGAAGCTAAATTAAGCTAAACAAGAG
CATGAAGGAAGACACAAGCCAGGAAGTACTGAAACCAGAGCTCTAGAAGGAAAGATGGGATGCTTACTGA
AGTTTTCAGGTGACTTGGATGACCAGACCTGTAGAGAAGATTTACTTCTCTCAAATCATGGTGA
AATAAAATGGGTGCACTTTGCCAGAGGACAAAAGAGGAATAATTCTCTCAAAGAAAAGGCTAAGGAA
GCACTTGAGAAAGCCAGAAATGCAAATAATGGTAACCTACTGTTAAGGAACAAAAGGTGACTTGGAAAG
TACTAGAAGGACATGCGGAAAAAGAAGCATTGAAAAAAATCACAGATGATCAGCAAGAATCCCTAAACAA
ATGGAAGTCAAAGGAGGCATGCAAGGTGAAGATTTAAAGGAAGCCATGTTTTTACAGCAGCCCGCAGA
TTTAAAGGAAAAGGAAGGGTAATAGGCCTGGTTATGCTGGGACCCAAAGGAAGGACAGTTTCATG
GAAGGAGGACAAGATTTGATGATGATCGTCGTCGTTGACCAATGAAAAGAGGAAGAGATGGAAGAGA
CAGAGAAGAACC CGCATCAAAGCATAAGAAAAGAGAAAATGGTCTAGAGACAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >MG206549 representing NM_009278
 Red=Cloning site Green=Tags(s)

MAENGDNKMTALEAKICHQIEYYFGDFNLPRDKFLKEQIKLDEGWVPLETMIKFNRLNRLTTDFNVIVQ
 ALSKSKAKLMEVSADKTKIRRSRPLPEVTDEYKNDVKNRSVYIKGFPTDATLDDIKEWLDDKGQILNI
 QMRRTLHKTFKGSIFAVFDSIQSAKKFVEIPGQKYKDTNLLILFKEDYFAKKNEERKQSKVEAKLKAKQE
 HEGRHKPGSTETRALEGKMGCLLKFSGDLDDQTCREDLHFLFSNHGEIKWVDFARGAKEGIIILFKEKAKE
 ALEKARNANNGNLLLRNKKVTWKVLEGHAEKALKKITDDQQESLNKWKSKGGHAGGRFKGSHVFTAARR
 FKGGKGNRPGYAGAPKGRGQFHGRRTRFDDDDRRRGPMPKRGRDGRDREEPASKHKKRENGARDK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009278


ORF Size: 2085 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_009278.4, NP_033304.1</u>
RefSeq Size:	1961 bp
RefSeq ORF:	1248 bp
Locus ID:	20823
UniProt ID:	<u>P32067</u>
Cytogenetics:	2 40.95 cM
Gene Summary:	Binds to the 3' poly(U) terminus of nascent RNA polymerase III transcripts, protecting them from exonuclease digestion and facilitating their folding and maturation.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG206549