

Product datasheet for **MG206715**

Serpina3a (BC063325) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Serpina3a (BC063325) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Serpina3a
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206715 representing BC063325 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCTTCATTGCAGCTCTGGGGCTCCTGATGGTGGGGATCTGTCTGCTGTACCTACTGGGCAACTG
 CAGATGGGCAACTGGGAGGGCACACTGCAGTCCAGAAAGACAGAGACCATGAAATACAACCTGGACAGTGT
 CACACTGGCTTCCATCAACACTGACTTTGCCTTCAGCCTCTACAAGAAGCTGGCTTTGGAGAATCCACAT
 AAAAATATTGTCTTCTCCCACTTAGCATCTCAGCTGCCTTGGCCCTCATGTCCCTGGGAGCAAAGGGCA
 ACACCTTGAAGAGATTCTAGAAGGTCTCAAGTTCAATCTCCAGAGACCCCTGAGGCAGACATCCACCA
 GAACTTTGGGCACCTCCTACAAATGCTCATTACGCCAGAGAACCAGGTTTCAAGTCAATGCAGGCAACGCC
 CTGTTTATTGATAAGCACCTGCAGATCTTGACAGAGTTCAAGGAGAAGGCAAGGGCTCTGTACAAGGCTG
 AGGCCTTCACAACAGACTTCCAGCGGCTCGTGAGGCCACAAAGCTCATCAATGACTATGTGAGGAAACA
 GACCCAGGGGAAGATCAAGGAAGTGTCTCAGACCTGCATAGGAACACATCCATGGCACTAGTGAATTC
 CTAACCTTTCAAGGCTTTTGAACGTGACTTTGACCCTGAAGACACGTTCTGGGAAATTTACATTGG
 ATAGGAAGAGGACAGTGAACGTGCCATGATGAAAACCGAGGAGCTAACCACAACTACTCCGGGACGA
 GGAGATGCAAAGCACTGTGGTGGAGCTGAATTACATAGGCAATGCCAGCTTCTGTTTCATCTCCCTGAC
 CAGGGCAGGATACAGCACGTGGAAGATAGCTTGAACCCAGTCTCTGAGGAAATGGAGGAAATCTTTGA
 GACCCAGGATGCTAGATGAACTCTCCCTGCCCAAGTTCTCCTTATCACAAGACTACAACCTAAATGACAT
 CCTTCCAGAGCTGGGTATCAAGGAAGTCTTCTCCACACAGGCTGACCTGTCTGGGATCACTGGAGCCAAA
 AATATAAGAGTCTCTCAGATGATCCACCAGGCTGCGCTGGATGTTACCGAGACACACAGAAGCAGATG
 TCATCACAATAGCCCGATATAACTTCCAATCTGCAAAAATTAAGGCTAAAATTGTTAAGGTGGACAGGGA
 ATTCTTGACCTTATTTAGACCAATGTTTAAGAGTATTAGCGTTATGGGTAAGGTTATCAACCTTTA
 ACAAAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >MG206715 representing BC063325
 Red=Cloning site Green=Tags(s)

MAFIAALGLLMVGICPAVTYWATADGQLGGHTAVQKDRDHEIQLDVTLASINTDFAFSLYKKLALENPH
 KNIVFSPLSISAALALMSLGAKGNTLEEILEGLKFNLPETPEADIHQNFHLLQMLIQPENQVQINAGNA
 LFIDKHLQILTEFKEKARALYKAEAFITDFQRPREATKLINDYVRKQTQGGIKELVSDLHRNTSMALVNF
 LNFQGGFWNVTFDPEDTFLGNFTLDRKRTVNVPMKTEELTTNYFRDEEMQSTVVELNYIGNASFLFILPD
 QGRIQHVEDSLQPQSLRKWRKSLRPRMLDEL.SLPKFSLSQDYNLNDILPELGIKEVFSTQADL.SGITGAK
 NIRVSQMIHQALDVTETHTEADVITIARYNFQSAKIKAKIVKVDREFLYLILDPMFKSISVMGKVINPL
 TN

TRTRPLE – GFP Tag – V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC063325

ORF Size: 1266 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [BC063325.1](#)


RefSeq Size: 1994 bp

RefSeq ORF: 1268 bp

Locus ID: 74069

Cytogenetics: 12 E

Product images:



Circular map for MG206715