

Product datasheet for **MG206727**

Cops3 (NM_011991) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cops3 (NM_011991) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cops3
Synonyms:	Csn3; Sgn3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206727 representing NM_011991 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGTCTGCTCTGGAGCAGTTCGTGAACAGTGTCCGACAGCTCTCAGCTCAAGGGCAAATGACTCAGC
TTTGTGAACATCAACAAGAGTGGAGAAGCTCCTTGCGAAGAAGTGTCCCATCTGGACTGTGCTTGG
GGCTTTAGATGTTCAAGAAGACTCTTTGGGTGCCTTGTCTGTTTGTGAAGTTCTCAATGCCAGT
GTTCTGACTTCGAAACACTATTCTCACAGTTCAGCTTTTCATCAGTACATGTAATGGCGAGCACATTC
GATATGCAACAGATACGTTTGTGGGCTTGGCCATCAGCTAACAAATGCCTTGTGAAAGAAAACAGCC
CCTTCGAGGAATTGGCATCCTTAAGCAAGCCATAGACAAGATGCAGATGAATACCAACCAACTGACCTCA
GTACATGCCGATCTCTGCCAGCTTTGTTACTAGCGAAGTGCTTTAAACCTGCCCTTCCATATCTTGATG
TGGATATGATGGATATCTGTAAGAGAATGGAGCCTATGATGCAAAACACTTTCTGTGCTACTACTATTA
TGGAGGAATGATCTATACGGGGCTGAAGAACTTTGAAAGAGCGCTGTACTTTTATGAGCAGGCTATAACT
ACTCTGCCATGGCCGTCAGCCATATCATGTTGGAATCATATAAAAAGTATTTCTAGTGTCTTTGATAT
TACTTGGCAAAGTACAACAGTGCCTAAATACACCTCTCAAATAGTGGGTAGATTTATTAAGCCTCTCAG
CAATGCTTACCATGAGTTAGCACAAAGTTTATTCACCAACAACCCGCTCTGAGCTGCGGAACCTGGTGAGC
AAGCACAGCGAGACCTTCACTCGAGACAACAACATGGGCCTGGTGAAGCAGTGCCTGTCTCTTTATA
AGAAGAACATTCAGAGGCTGACAAAGACTTTTTAACACTATCATTACAAGATATGGCAAGCCGTGTACA
GTTATCTGGACCACAGGAGGCGAGAAAAATATGTTCTACACATGATAGAAGATGGAGAAAATTTTTGAAGT
ATTAACGAGAAGGATGGAATGGTCAGTTTCCATGATAACCCTGAGAAGTATAATAACCCAGCCATGCTTC
ATAACATCGACCAAGAGATGCTAAAATGCATAGAGCTGGATGAGCGACTGAAGGCTATGGACCAGGAGAT
CACGGTGAACCCCGAGTTTGTCAAAGAGCATGGGCTCACAGAAGATGACTCAGGAAACAAGCCCTCC
AGCTACTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG206727 representing NM_011991
 Red=Cloning site Green=Tags(s)

MASALEQFVNSVRQLSAQGQMTQLCELINKSGELLAKNLSHLDTVLGALDVQEHSLGVLAVL FVKFSMPSP
 VPDFETLFSQVQLFISTCNGEHIRYATDTFAGLCHQLTNALVERKQPLRGIGILKQAIKMQMNTNQLTS
 VHADLCQLCLLAKCFKALPYLDVDMMDICKENGAYDAKHFLCYYYYGGMIYTGKKNFERALYFYEQAIT
 TPAMAVSHIMLESYKKYILVSLILLGKVQQLPKYTSQIVGRFIKPLSNAYHELAQVYSTNNPSELRLNLS
 KHSETFTRDNNMGLVKQCLSSLYKKNIQRLTKFTLTLSLQDMASRVQLSGPQEAKEYVLHMIEDGEIFAS
 INQKQGMVSVFHDNPEKYNNPAMLHNIDQEMLKCIELDERLKAMDQEITVNPQFVQKSMGSQEDDSGNKPS
 SYS

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_011991

ORF Size:	1269 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_011991.1 , NP_036121.1
RefSeq Size:	1583 bp
RefSeq ORF:	1272 bp
Locus ID:	26572
UniProt ID:	O88543
Cytogenetics:	11 B1.3
Gene Summary:	Component of the COP9 signalosome complex (CSN), a complex involved in various cellular and developmental processes. The CSN complex is an essential regulator of the ubiquitin (Ubl) conjugation pathway by mediating the deneddylation of the cullin subunits of SCF-type E3 ligase complexes, leading to decrease the Ubl ligase activity of SCF-type complexes such as SCF, CSA or DDB2. The complex is also involved in phosphorylation of p53/TP53, c-jun/JUN, IκappaBα/NFKBIA, ITPK1 and IRF8/ICSBP, possibly via its association with CK2 and PKD kinases. CSN-dependent phosphorylation of TP53 and JUN promotes and protects degradation by the Ubl system, respectively. Essential to maintain the survival of epiblast cells and thus the development of the postimplantation embryo.[UniProtKB/Swiss-Prot Function]