

Product datasheet for **MG206888**

Mrps5 (NM_029963) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrps5 (NM_029963) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Mrps5
Synonyms:	1620401I16Rik; AI850294
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG206888 representing NM_029963 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCGGCCGTACGCGCTGCAGGCTGTCTCCCTGCATTGTGTAGCTTGCAGGCAGGTCATTTCTGT
CTAGGCAACTCTCTTTAAACGCCTTTCCAGTAGCAGCCACCTCCTTTTGGCAGTGAAGACAGCTCTTAG
CCACGGCTCATTGTCATCCAGGGAACAAGGCGCAACCACTGTCTGACCAGCTTGAGCCATGTGCTACAG
ACACAGTGCTGCGTTTCTCTCCCGCACTGGACGGGCCAGCAGTGCCGGCCCTACAGCTTCTTACCA
AACTGACTGCAGAAGAGCTCTGAAAGGCGCATTAGCAGAGACTGGAGCTGGAGCAAGAAAGGCAGAGG
CAAAAGAACAAAGAAAAAGAGGAAGGATTTGAACAGGGGCCAGATCATTGGTGAAGGGCGCTCTGGC
TTCTGTGGCCTGGTTTGAACGTTCTCTGATAAAAAAGTGGGGTCGTCCAGAATCGGCCAGAGAAGCA
AAGAGGAGCAGCAGAAGGTGGAGGCCACCATGGTCGAGCAGCGAGAGGAGTGGACCGGAAGAGGAAGAT
CAAAGTTAAAGGGAGCGCGGCTGGAGTGGGAACAGTGGGGGGTGTGAGTATTGGCCCCCAGACCCG
GGACCAATGGAGAAACATATGAAGACTTTGATACCAGGATTCTTGAGGTAAGAAATGTCTTCAATATGA
CAGCAAAAGAGGGAAGAAAGAGTCGGTCCGAGTCCTGGTTGCTGTTGGGAATGGAAATGGAGCTGCAGG
TTTTGCTATTGGGAAAGCTGTGACCGGGGAGACGCTTTCAGAAAGGCAAGAACCGAGCGATTCTTAT
TTGCATTACATAGAACGATACGAAGGACATACAATATTCATGACATTTCCCTAAGATTCAAAGGACAC
AGATCCGGATGAAGAAACAACCCAGAGGTTACGGCCTCCGCTGCCACAGGGCCATTATCACCATCTGCCG
CCTCATCGGCATCAAGGACATGTATGCGAGGGTCACCGGGTCCATGAACATGCTCAACCTCACCCGGGGC
CTCTTCCATGGGCTTGCCCGCCAGGAAACCCACCAACATCTGGCGGATAAGAAGGGTCTCCACGTGGTGG
AATTCCGGGAGGAATGTGGGCTCTGCCTATCGTGGTGGCCTCCACATGGGGCCTTGAAGAGGCC
GGAGCCGAACCTGAGGTTCTGTATACCAAGCTGGACTGGCAGGACGTGAAGGCCATGCAGGACTGAAG
CGCTCTGTGTGGTTAATTTAAAGAGGCTGCCACC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG206888 representing NM_029963
 Red=Cloning site Green=Tags(s)


MAAAVRAAGCLPALCSLQAGHFLSRQLSLNAPFVAATSFLAVKTALSHGSLSSRETRRNHCLTSLSHVLQ
 TQCCVSSPGNWTGQCRPYSFFTKLTAEELWKGALAEAGARKGRGKRTKKKKRDLNRGQIIGEGRS
 FLWPGLNVPLIKSGVVQNIQGRSKEEQKVEATMVEQREEWDRKRKIKVKRERWSGNTWGGVSI
 GPPDP
 GPNGETYEDFDTRILEVRNVFNMTAKEGRKKSVRVLVAVGNGNGAAGFAIGKAADRGDAFRKAKNRAIHY
 LHYIERYEGHTIFHDISLRFKRTQIRMKKQPRGYGLRCHRAIITICRLIGIKDMYARVTGSMNMLNLTRG
 LFHGLARQETHQLADKKGLHVVEFREEGGLPIVVASPHGALSKEPEPEPEVPTKLDWQDVKAMQGLK
 RSVWFNLKRPAT

TRTRPLE – GFP Tag – V

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_029963

ORF Size: 1296 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_029963.2](#), [NP_084239.1](#)

RefSeq Size: 1448 bp


RefSeq ORF: 1299 bp

Locus ID: 77721

UniProt ID: [Q99N87](#)

Cytogenetics: 2 F1

Product images:



Circular map for MG206888