

Product datasheet for **MG207287**

Rnh1 (NM_145135) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rnh1 (NM_145135) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Rnh1
Synonyms:	AW546468; C80305; RNH
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG207287 representing NM_145135
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGTCTTGACATCCAGTGTGAGCAGCTGAGTGTGCCCGTGGACAGAGCTCCTTCTCTGATCCAGC
 AATACGAAGTGGTCAGGCTGGATGACTGTGGCCTCACTGAAGTGGGTGCAAAGACATCAGCTCAGCAGT
 CCAAGCCAACCCTGCCCTGACAGAGCTCAGCCTTCGCACCAATGAACTGGGTGATGGCGGTGTGGGTCTG
 GTGCTCCAGGGCCTGCAGAATCCCACTTGAAGATCCAGAAGCTGAGCCTTCAGAACTGCGGCTTACCGG
 AGGCTGGCTGTGGGATCCTGCCTGGCATGCTGCGCTCTTGTCTACCCTGCGTGAGCTACATCTCAATGA
 TAACCCTATGGGGGACGCAGGCCTGAAGCTGCTCTGTGAAGGACTTCAGGACCCCAAGTCCGCTTGTGAG
 AAGCTTCAGTTGGAATACTGTAACCTCACAGCTACCAGCTGCGAGCCCTGGCTTCAGTGTCTCAGGGTGA
 AAGCTGACTTTAAGGAGTAGTATTGAGCAACAACGACCTCCATGAACCTGGTGTCCGTATTCTGTGCCA
 GGGCCTGAAGGACTCTGCCTGTCAACTGGAGTCACTCAAGCTGGAGAAGTGTGGTATCACAGCAGCCAAC
 TGCAAGGATCTGTGCGATGTTGTGGCCTCAAAGCTTCGCTACAGGAACTGGACTTGAGCAGCAACAAGC
 TAGGTAATGCAGGCATTGCAGCACTGTGCCAGGACTGTGCTTCCCAGCTGTAAGCTCAGGACTCTGTG
 GCTCTGGGAATGTGATATCACTGCAGAGGGCTGCAAGGACCTGTGCCGTGTCTCAGAGCCAAGCAGAGC
 CTGAAGGAACTCAGCCTTGCTTCCAATGAGCTGAAGGATGAGGGTGCCCGGCTGTGTGTGAGAGCCTAT
 TAGAGCCTGGCTGTGAGTGGAGTCACTGTGGATAAAGACCTGTAGCCTCACAGCTGCCTCCTGTCCCTA
 CTTCTGCTCCGTGTTGACAAAAGCCGTTCTCTGTTGGAGCTGCAATGAGCAGCAACCCACTGGGAGAC
 GAGGGAGTCCAGGAGCTTTGCAAGGCCCTGAGCCAGCCAGATACAGTGTGCGTGAGCTTTGGTTGGGAG
 ACTGTGATGTGACAAATAGTGGTTGCAGCAGCCTTGCCAATGTCCTGCTGGCCAACCGCAGCCTGAGGGA
 ACTGGACCTCAGTAACTGCACTGGGGGCCAGGTGTCTGCAACTGCTGGAGAGCCTCAAACAGCC
 AGCTGCACCCCTCAGCAGCTCGTCTGTATGACATTTACTGGACAAATGAGGTGGAAGAAGCAGCTTCGGG
 CCCTGGAGGAGGAAAGCCCTTCCCTGAGGATCATTTCC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG207287 representing NM_145135
 Red=Cloning site Green=Tags(s)

MSLDIQCEQLSDARWTELLPLIQYEVVRLDDCGLTEVRCKDISSAVQANPALTELSLRTNELGDGGVGL
 VLQGLQNPTCKIQKLSLQNCGLTEAGCGILPGMLRSLSTLRELHLNDNPMGDAGLKLCEGLQDPQCRL
 KLQLEYCNLTATSCPELASVLRVKADFKELVLSNNDLHEPGVRIILCQGLKDSACQLESKLENCGITAAN
 CKDLCDVVASKASLQELDLSSNKLGNAGIAALCPGLLLPCKLRTLWLWECIDITAECKDLCRVLRKQS
 LKELSLASNELKDEGARLLCESLLEPGCQLESLEWIKTCSLTAASCPYFCSVLTKSRSLLELQMSNPLGD
 EGVQELCKALSQPDTVLRRELWLGDCDVTNSGSSLANVLLANRSLRELDLSNNCMGGPGVQLLESLKQP
 SCTLQQLVLYDIYWTNEVEEQLRALEEERPSLRIS

TRTRPLE – GFP Tag – V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_145135

ORF Size: 1368 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145135.3](#), [NP_660117.2](#)

RefSeq Size: 1698 bp

RefSeq ORF: 1371 bp

Locus ID: 107702

UniProt ID: [Q91VI7](#)

Cytogenetics: 7 F5

Gene Summary: Ribonuclease inhibitor which inhibits RNASE1, RNASE2 and ANG. May play a role in redox homeostasis (By similarity).[UniProtKB/Swiss-Prot Function]