

Product datasheet for **MG207364**

Trim61 (BC064738) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trim61 (BC064738) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Trim61
Synonyms:	2czf61; AU018704; C85201; E330039K03Rik; Rnf35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG207364 representing BC064738
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGCTTGCTGCAGGGCTGGCCGATCTCCAAGCAGAGGCTCGTGTTCTGTCTGTAAGGAATACCTCA
 AGAACCCAGGGACCATGGAGTGTGGCCATAAAGTCTGTCTGTCTGCACCAAGTGTGTTCTGGGAGGATCT
 AAAGGGTAGTGTCCCCTGCCCTCTTGTCATTTGAACTGTCCAGGAAGGAACCTTTTGTAGTAATCAGCCA
 CCAGGTGAGTTCACCGAAGTTGCCATCTTGCCCTGAAAGAACTAGCAGGAGTGAAGGCAAGAAGGAG
 AGCATCCATGGGAAAAGCCACAACAGGTTTCTGTCCACTGCGCAGCTTCCCCCGAATCACTTTGTTTG
 GCCCGTGGAAAATGCTATCGTCTGCCACAGGAAACAAATTGACTGTTACATTAAGCTGTGGAGGGAGATG
 GTGGAGCCAATCCAGAAAGCAATTAGCACACAAAGGAAAGGATCGCTAGAACTGAGAAAAAGGACAGAAT
 GCAGAAGAGAAAGAGTGAGGTCTCAATATGAACAGAATAGGTTGTTCTCCAGGTTGAACGGAAGAAGAT
 TCTTGGGCACTTACAAAGTGAAGAGATGGATGGGCTATCCAACTAAACAAAAACCTAGCGAAATTATCA
 GATCACACCTCCTTGTTAAAGAGTCTGTTGAAAGATCTAAAGAGCCAGTGTGAAAAGTCAGAGCTGGCAT
 TGCTGGCAAGTGTTAAAGATTGGTACGAGAGATACAAACAGTCGGGATGGCCTGAGATGGCCTTAGTAAC
 ATTGAAAGGTACGACTACAGGTTACCTCTACAGTATTCTGACTTAGACAGAATTATCAAGCGATTTCAA
 GTTGATGTCTTCTGGATCCCGAAACAGCAAACTATAAGCTTACAGTCTCAAAGGATAAAAAACAGTTA
 AATATGGGAGTTGGGAGAGGGTGCCCTATAGCTCAAGAAAATTTGTTTTGACCTGCTGTTCTGGGATC
 TGAGGGATATAGTTCTGACAGACAATACTGGGAGGTTGATGTGGAACCAAGATGGAATGGGTTCTAGGG
 GTCTGTAGGGAGCCTTTTCCAGGAGCAGGTATGCCGAACGGTTCTCTGATGAACAGTTCTCAGTGCAGG
 ATGGATTATGGGAGTTGGGCTAGACAATCCCAAACCTTATGTTGCATTGGGTGGTGAGAAAAATCCATCT
 TCTTCCCAAAGTATCACCTACTAGGATTGGGATTTTTTTAGATTCTGAAATGTATGAGGTTTCGTTTTAC
 AATTTGAGGGATAAATCCCTGCTGTACAGGTTTAGCGACTTACCTCGGGGCACATTTTGGCCTTACTTCT
 ACACTGGATATGACTACGAGCCTCTGAAGATTGACAGTTGCAGATCTTGAG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:


>MG207364 representing BC064738
 Red=Cloning site Green=Tags(s)

MELAAGLADLQAEARCSVCKEYLKNPGTMECGHKVCLSCTSVFWEDLKGSVPCPSCHLNCPGRNFLSNQP
 PGFEFEVAILPPERTSRSGRQGEHPWEKQQVSCPLRSFPNHFVWPVENAIVCHRKQIDCYIKLWREM
 VEPIQKAISTQRKGSLELRKTECRRERVRSQYEQNRLFQVERKKILGHLQSEEMDGLSKLNKNLAKLS
 DHTSLLKSLKDLKSQCEKSELALLASVKDWYERYKQSGWPEMALVTLKGYDYRLPLQYSDLDRIKRFQ
 VDVILDPETANYKLTVSKDKKTVKYGSWERVYPYSSRKFCFDPVGLGSEGYSSDRQYWEVDVELKMEWVLG
 VCREPFRSRYAERFSDEQFSVQDGLWGVGLDNSQTYVALGGEKIHLPLPKVSPTRIGIFLDSEMYEVSFY
 NLRDKSLLYRFSDLPRGTFWPYFYTYDYEPLKICTVADLE

TRTRPLE – GFP Tag – V

Restriction Sites:


SgfI-MluI

Cloning Scheme:


ACCN:	BC064738
ORF Size:	1383 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC064738 , AAH64738
RefSeq Size:	2107 bp
RefSeq ORF:	1385 bp
Locus ID:	260296

Cytogenetics: 8 B3.1

Product images:



Circular map for MG207364