

Product datasheet for **MG207626**

Gdf10 (NM_145741) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gdf10 (NM_145741) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Gdf10
Synonyms:	Bmp3; Bmp3b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG207626 representing NM_145741
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCTCCAGGTCTGCTCGGATCAGCTTGGGGTCCCAGCTGCTGCCATGGTGCCGCTGCTCCTGCTGC
 TCGGGGGCGCAGGCTGCGGCCACAGGGGCCCTCATGGTCTCATTGCCCTCGGCAGCTGCCGGTCTGCA
 GGGGACAGGGACTCCCAGCAGTCAACCGGGGACGACAGCCGCTCTGGGCCAGGCGCCAGGACATG
 GTCGCTATCCACATGCTCAGGCTCTATGAGAAGTACAACCGAAGAGGTGCTCCACCGGGAGGAGCAACA
 CCGTCCGAAGCTTCCGTGCCCGGCTGAAATGATCGACAAAAGCCTGTGTATTTCTTCACTTGACTTC
 CATGCAAGACTCAGAAATGATCCTCACAGCCGCTTCCACTTCTACTCAGAACCTCCACGGTGGCCCGG
 GCTCGTGAGGTATTCTGCAAGCCCCGAGCTAAGAACGCATCCTGCCGCCTCTGACCCAGGGCTGCCTG
 CACGCTTGACCTAATCTCCGCAGTCTTCCAGAACCCGCCACTCAGGGGCTGCTCCGCGGGGCCAT
 GGCCTGACGCCTCCACCAGTGGCCTGTGGCAGGCCAAGGACATCTCCTCAATCATCAAGCTGCCGA
 AGGGATGGAGAGCTGTTCTCTGCTCAGCTGGATACTGGGGAGAAGGACCCCGGAGTGCCACGGCCCA
 GTTCCACATGCCCTATATCCTTGTCTACGCCAATGACCTGGCCATCTCCGAACCAACAGTGTAGCAGT
 GTCGCTACAGAGATACGACCCATTTCCAGCTGGAGACTTTGAGCCTGGAGCAGCCCCAACAGCTCAGCT
 GATCCCCGCGTGCAGGGCGGCTCAGGTGTCAAACCCCTGCAAGACAATGAATGCCGGGGCTGGATG
 AAAGACCAGCGCTGCCCTGCATGCCAGAAATTTCCACAAGCAGCAGTTCCTGGTCCAGTCTTTCCGGGC
 ACTGAAACCCCGCAGGGGGCGAAAGACCGCAAGAAGAAGGACCAGGACACATTCACCGCCGCTCCTCT
 CAGGTGCTGGACTTTGACGAGAAGCAGATGCAGAAAGCCAGGAGCGGCAGTGGGATGAGCCCCGGTCT
 GCTCCAGGAGTACCTGAAGGTGATTTTGCAGACATCGGGTGGAAATGAATGGATCATCTCTCCAAATC
 CTTTGACGCCTACTACTGTGCTGGGCTGCGAGTTCCTCCATGCCCAAGATTGTCGCCCATCCAACCAT
 GCCACCATCCAGAGCATCGTCAGAGCTGTGGCATTGTCCTGGCATCCCAGAGCCATGCTGTGTTCCAG
 ACAAGATGAACTCCCTTGAGTCTTTTCTGGATGAAAATCGGAATGCGGTTCTGAAGGTGTACCCCAA
 TATGTCCGTAGAGACCTGTGCTGTCGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG207626 representing NM_145741
 Red=Cloning site Green=Tags(s)

MAPGPARISLGSQLLPMVPLLLLLRGAGCGHRGPSWSSLPSAAAGLQGDRDSQQSPGDAAAALGPGAQDM
 VAIHMLRLYEKYNRRGAPPGGGNTVRSFRARLEMIDQKPVYFFNL TSMQDSEMIL TAAFHFYSEPPRWPR
 AREVFCKPRAKNASCRLLTPGLPARLHLIFRSLSQNTATQGLLRGAMALTPPPRGLWQAKDISSIIKAAR
 RDGELLLSAQLDTGEKDPGVPRPSSHMPYILVYANDLAISEPNSVAVSLQRYDPPFAGDFEPGAAPNSSA
 DPRVRAAQVSKPLQDNELPGLDERPAPALHAQNFHKHEFWSSPFRALKPRTGRKDRKKKQDFTTAASS
 QVLD FDEKTMQKARRRQWDEPRVCSRRYLKVDFA DIGWNEWII SPKSFDAYYCAGACEFMPKIVRPSNH
 ATIQSIVRAVGIVPGIPEPCVPDKMNSLGVLF LDENRNAV LKVYPNMSVETCACR

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	<p>This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.</p>
Components:	<p>The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).</p>
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145741.3
RefSeq Size:	2613 bp
RefSeq ORF:	1431 bp
Locus ID:	14560
UniProt ID:	P97737
Cytogenetics:	14 20.8 cM
Gene Summary:	<p>This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. This protein has been shown to promote neural repair after stroke and may act as a tumor suppressor. [provided by RefSeq, Jul 2016]</p>