

Product datasheet for **MG208087**

Hsf1 (NM_008296) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hsf1 (NM_008296) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Hsf1
Synonyms:	AA960185; Hsf1alpha; Hsf1beta; HSTF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG208087 representing NM_008296
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGATCTGGCCGTGGGCCCGGTGCAGCGGGGCCAGCAACGTCCCGGCCTTCTAACCAAGCTGTGGA
 CCCTCGTGAGCGACCCGACACAGACGCGCTCATCTGCTGGAGCCGAGTGGGAACAGCTCCACAGTGT
 TGACCAGGGCCAGTTTGCCAAGGAGGTGCTGCCAAGTACTTCAAGCACAAACATGGCTAGCTTCGTG
 CGGCAGCTCAACATGTATGGCTCCGAAAAGTAGTCCACATTGAGCAGGGTGGCCTGGTCAAGCCTGAGA
 GAGATGACACCGAGTTCAGCATCCTTGTCTTCTGCGTGGACAGGAACAGCTCCTTGAGAACATCAAGAG
 GAAAGTGACCAGCGTGTCCACCCTGAAGAGTGAGGACATAAAAAACGCCAGGACAGTGTACCCCGGCTG
 TTGACAGATGTGCAGCTGATGAAGGGGAAACAGGAGTGTATGGACTCCAAGCTCCTGGCCATGAAGCAG
 AGAACGAGGCCCTGTGGCGGGAGGTGGCCAGCCTTCGGCAGAAGCATGCCAGCAGCAAAAAGTTGTCAA
 CAAGCTCATTAGTTCCTGATCTCACTGGTGCAGTGAACCGGATCCTGGGGTGAAGAGAAAGATCCCT
 CTGATGTTGAGTGACAGCAACTCAGCACACTCTGTGCCAAGTATGGTCCGACAGTACTCCCTGGAGCATG
 TCCATGGTCTGGCCATACTCAGCTCCATCTCCAGCCTACAGCAGCTCTAGCCTTTACTCCTCTGATGC
 TGTACCAGCTCTGGACCATAATCTCCGATATCACTGAGCTGGCTCCCACCAGCCCTTTGGCCTCCCCA
 GGCAGGAGCATAGATGAGAGGCCCTGTGCCAGCAGCACTCTGGTCCGTGTCAAGCAAGAGCCCCCAGCC
 CACCTCACAGCCCTCGGGTACTGGAGGGCAGCCCTGGCGGCCATCCTCCATGGATACCCCTTTGTCCCC
 AACTGCCTTATTGACTCCATCCTTCGAGAGAGCGAGCCTACCCCTGTGCCTCAAACACAGCCCTATG
 GACACAACCGGAGCCCAAGCCCCGCACTCCCAGCCCTCCACCCTGAGAAGTGCCTCAGCGTAGCCT
 GCCTAGACAAGAACGAGCTAAGTATCACCTGGATGCCATGGACTCCAACCTGGACAACCTGCAGCAT
 GCTGACAAGCCACGGCTTCAAGTGTGGACACCAGTGCCTGCTGGACATTACAGGAGCTTCTGTCTCCAAA
 GAGCCTCCAGGCCTATTGAGGCAGAGAACAGTAACCCCGACTCAGGAAAGCAGCTGGTGCCTACACGG
 CTCAGCCTCTGTTCTGCTGGATCCTGATGCTGTGGACACAGGGAGCAGTGAAGTGCCTGTCTTTGA
 GCTGGGGGAGAGCTCCTACTTCTCTGAGGGGATGACTACACGGATGATCCCACCATCTCTCTTCTGACA
 GGCCTGAACCCATAAAGCCAAGGACCCCACTGTCTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:


>MG208087 representing NM_008296
 Red=Cloning site Green=Tags(s)

MDLAVGPGAAGPSNVPAFLTKLWTLVSDPDTDALICWSPSGNSFHVFDQGFQFAKEVLPKYFKHNNMASFV
 RQLNMYGFRKVVHIEQGLVKPERDDTEFQHPCFLRGQEQLENIKRKVTSVSTLKSEDIKIRQDSVTRL
 LTDVQLMKGKQECMDSKLLAMKHENEALWREVASLRQKHAQQQKVVNKL IQFLISLVQSNRILGVKRIKIP
 LMLSDSNSAHSVPKYGRQYSLEHVHGPYPYSPAYSSSLYSSDAVTS SGP IISDITELAPTSPLASP
 GRSIDERPLSSSTLVRVKQEPSPPHSPRVLEASPRPSSMDTPLSPTAFIDSILRESEPTAASNTAPM
 DTTGAQAPALPTPSTPEKCLSVACL DKNELSDHLDAMDSNLDNLQTMLTSHGFSVDTSALLDIQELLSPQ
 EPPRPIEAENSNPDSGKQLVHYTAQPLFLDDPAVDTGSSELVLFELGESSYFSEGDDYTDPTISLLT
 GTEPHKAKDPTVS

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Plasmid Map:


ACCN: NM_008296

ORF Size: 1509 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008296.3](#)

RefSeq Size: 2018 bp

RefSeq ORF: 1512 bp

Locus ID: 15499

UniProt ID: [P38532](#)

Cytogenetics: 15 35.95 cM

Gene Summary: Function as a stress-inducible and DNA-binding transcription factor that plays a central role in the transcriptional activation of the heat shock response (HSR), leading to the expression of a large class of molecular chaperones heat shock proteins (HSPs) that protect cells from cellular insults' damage. In unstressed cells, is present in a HSP90-containing multichaperone complex that maintains it in a non-DNA-binding inactivated monomeric form. Upon exposure to heat and other stress stimuli, undergoes homotrimerization and activates HSP gene transcription through binding to site-specific heat shock elements (HSEs) present in the promoter regions of HSP genes. Activation is reversible, and during the attenuation and recovery phase period of the HSR, returns to its unactivated form. Binds to inverted 5'-NGAAN-3' pentamer DNA sequences. Binds to chromatin at heat shock gene promoters. Plays also several other functions independently of its transcriptional activity. Involved in the repression of Ras-induced transcriptional activation of the c-fos gene in heat-stressed cells. Positively regulates pre-mRNA 3'-end processing and polyadenylation of HSP70 mRNA upon heat-stressed cells in a symplekin (SYMPK)-dependent manner. Plays a role in nuclear export of stress-induced HSP70 mRNA. Plays a role in the regulation of mitotic progression. Plays also a role as a negative regulator of non-homologous end joining (NHEJ) repair activity in a DNA damage-dependent manner. Involved in stress-induced cancer cell proliferation in a IER5-dependent manner.[UniProtKB/Swiss-Prot Function]