

Product datasheet for **MG208452**

Cyp4f15 (BC021377) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyp4f15 (BC021377) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cyp4f15
Synonyms:	AI787289; Cyp4f-15
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG208452 representing BC021377
 Red=Cloning site Blue=ORF Green=Tags(s)
 TTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGGGTTTTTCAGGATGCCACAGTTGGACCTGTCTTGGCTGGGACTGAGGCTGGAGGCATCCTCACCAT
 GGCTGCTTTTGTGCTGATCGGAGCTTCTGGCTCTTGGCCCGTGTCTGACTCAGACCTATATCTTCTA
 CAGAACATATCATCATCTTTGTGATTTCCCTCAACCCCCAAATGGAATTGGTTCTTGGGTACCTAGGC
 ATGATCACACCCACTGAGCATGGTCTGAAGGAGGTGACTAACCTGGTGGCCACCTACCCCGAGGCTTCA
 TGACATGGCTGGGTCTATCATTCTATCATCACTCTGTGCCACCCTGACATCATTGCATCTGTTCTCAA
 TGCCTCAGCTTCGGTAGCCCTAAAGGAAGTAGTCTTCTACAGTTTCTGAAACCTTGGCTGGGGGATGGG
 CTCTTGCTGAGTGATGGTGACAAGTGGAGCAGCCACCGTCGCATGTTGACACCTGCTTTCCATTTCAACA
 TCCTAAAACCTATGTGAAGATTTCAATGACAGCACCAACATAATGCATGCCAAGTGGCAGCACCTGGC
 CTCAGGTGGCAGCGCCCTCTGGATGTGTTTGAAGACATCAGCCTCATGACCTTGGACAGCCTGCAGAAA
 TGTGTCTTCAGTTTTGATAGCAACTGCCAGGAGAATCCAGTGAGTATATTTCTGCCATCTTGGAGCTCA
 GTGCCCTAGTGACCAACGGTACCACCACTGCTTCTACACATAGACTCACTCTATCAGCTCACCTGTAG
 TGGGCGGCGCTTCCACAAGGCTGCCACCTGGTGACAGCTTACAGATGCTGTCATCCAGGATCGACGA
 CGTACTCTCCCGAGTAAGCATGAGGATGATGTCCTCAAGGCCAAGGCTAAGTCCAAAACCTCTGGATTTCA
 TTGATGTGCTGCTGCTGAGCAAGGATGAAGATGGAAGGAGCTGTGAGATGAGGACATTAGAGCAGAGGC
 TGACACCTTCATGTTTGAAGGCGCATGACACCACAGCCAGCGGGCTCTCCTGGATCCTGTACAATCTGGCA
 AGGCACCCAGAATACCAGGAGCGCTGCCGCGAGGAGGTGCAGGAGCTGCTGAGGGATCGAGAGTCTACAG
 AAATTGAGTGGGACGACTTGGCCAGCTACCTTCTGACCATGTGCATCAAGGAGAGTCTGCGGCTGACG
 TCCCCAGTCACGGTCATTTCCCGTCGTTGTACCCAAGACATTGTGCTCCAGATGGCCGGGTCAACCT
 AAAGGTGTATCTGCATCATCAATATTTTCGCGACCCATCACAACCCAACTGTGTGGCCTGACCCTGAGG
 TTTATGACCCCTTCCGCTTTGACCCAGAGAATATCAAGGACAGGTACCTCTGGCTTTTATTCCTTCTC
 CGCTGGGCCCAGGAAGTGCATAGGACAGACGTTTGCTATGAACGAGATGAAGGTGGCGCTGGCGTGACG
 CTGCTGCGCTTCCGTGTGCTGCCAGATGACAAGGAGCCGCGCGGAAGCCGAGTTGATCTGCGCGCAG
 AGGGCGGGCTGTGGCTGCGGGTGGAGCCGCTGAGCACACAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG208452 representing BC021377
 Red=Cloning site Green=Tags(s)
 MGFFRMPQLDLSWLGRLLEASSPWLLLLLIGASWLLARVLTQTYIFYRTYHHLCDFPQPPKWNWFLGHLG
 MITPTEHGLKEVTNLVATYPQGFMTWLGPIIPIITLCHPDIIIRSVLNASASVALKEVVFYSFLKPWLGDG
 LLLSDGDKWSSHRRMLTPAFHFNLKPYVKIFNDSTNIMHAKWQHLASGGSARLDVFENISLMTLDSLQK
 CVFSFDSNCQENPSEYISAIELLSALVTKRYHQLLLHIDSLYQLTCSGRRFHKACHLVHSFTDAVIQDRR
 RTLPSKHEDDLKAKAKSKTLDLIDVLLLSKDEDGKELSDIEDIRAEADTFMFEGHDTTASGLSWILYNLA
 RHPEYQERCRQEVQELLRDRESTEIEWDDLAQLPFLTMCIKESLRLHPPVTVISRRCTQDIVLPDGRVIP
 KGVICIIINIFATHHNPTVWPDPEVYDPFRFDPENIKDRSPLAFIPFSAGPRNCIGQTFAMNEMKVALALT
 LLRFRVLPDDKEPRRKPELILRAEGLWLRVEPLSTQ

TRTRPLE - GFP Tag - V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: BC021377

ORF Size: 1581 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC021377 , AAH21377
RefSeq Size:	2120 bp
RefSeq ORF:	1583 bp
Locus ID:	106648
Cytogenetics:	17 B1

Product images:



Circular map for MG208452