

Product datasheet for **MG208529**

Mgrn1 (NM_029657) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Mgrn1 (NM_029657) Mouse Tagged ORF Clone
 Tag: TurboGFP
 Symbol: Mgrn1
 Synonyms: 2610042J20Rik; md; mKIAA0544; nc
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-AC-GFP (PS100010)
 E. coli Selection: Ampicillin (100 ug/mL)
 Restriction Sites: SgfI-MluI
 Cloning Scheme:

Cloning sites used for ORF Shuttling:




ACCN: NM_029657
 ORF Size: 1596 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_029657.4
RefSeq Size:	3138 bp
RefSeq ORF:	1599 bp
Locus ID:	17237
UniProt ID:	Q9D074
Cytogenetics:	16 2.48 cM
Gene Summary:	E3 ubiquitin-protein ligase. Mediates TSG101 monoubiquitination at multiple sites. Plays a role in the regulation of endosome-to-lysosome trafficking. Impairs MC1R- and MC4R-signaling by competing with GNAS-binding to MCRs and inhibiting agonist-induced cAMP production. Does not inhibit ADRB2-signaling. Does not promote MC1R ubiquitination (By similarity). Acts also as a negative regulator of hedgehog signaling (PubMed:29290584). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG208529