

Product datasheet for **MG208728**

Slc9a7 (BC058750) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc9a7 (BC058750) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Slc9a7
Synonyms:	A530087D17Rik; NHE7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG208728 representing BC058750
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGAGCCTAGTGACGCGCGCGCCCGGGTCCCGGAGAGCTTTTCGGGGGCTGTCGCCGCGCTGCTGC
 TGCTGCCACTGTTGCCTGTGCTGCTGGGTGCGGGGCTGCGGGCTGGGGCCGCGGCCTCTCTGGGAGGC
 GGCTGAGGACAGCAGCGCCATGGAGGAGCTTGCCACCGAGAAGGAGGCGGAGGAGAGCCACCGGCAGGAC
 AGCGTGAGCCTGCTCACCTTCATCCTGCTGCTCACCTCACCATCTCAGCATCTGGCTTTCAAGCACC
 GCCGGGTGCGTTTCCTTCATGAGACCGGCTGGCCATGATTTATGGGCTCATTGTTGGAGTGATCCTGAG
 ATATGGCACCCCTGCTACCAGCGGTGATGACAAGTCACTCAGCTGTACTCAGGAAGACAGAGCCTTTAGC
 ACTTTATTAGTGAATGTCAGTGGGAAGTTCTTTGAATACACTCTGAAAGGAGAAATCAGCCCTGGCAAGA
 TCAACAACGTGGAGCAGAATGACATGCTGCGGAAGGTAACATTTGATCCAGAGGTATTTTCAACATTCT
 GCTGCCTCCAATTATTTCCATGCTGGATACAGCTTAAAGAAGAGACACTTTTTTCAGAAATCTTGGGTGCG
 ATTCTGGCCTATGCTTTCTTGGGACTGCTGTTTCTTGCTTCATTATTGGAATCTTATGTATGGAGTGG
 TGAAACTCATGAAGATTGTGGGCAGCTCTCTGATAAATTTTACTACACAGATTGCTCTTTTTTGGAGC
 AATTATTTCTGCCACTGACCCAGTGACTGTGCTGGCAATCTTTAACGAACACATGCAGATGTGGATCTT
 TATGCTCTTCTGTTTGGAGAGAGTGCTGAATGATGCAGTTGCCATTGTACTATCCTCGTCTATTGTTG
 CCTATCAACCAGCAGGACTAAACACTCATGCGTTTGATGCTGCTGCTTTTTTAAGTCAGTTGGCATTCT
 CCTAGGTATATTTAGTGGCTCTTTTACCATGGGAGCTGTAAGTGGTGTGTTACTGCTCTAGTGACCAAG
 TTTACTAAGCTGCACTGCTTCCCTGCTGGAGACAGCGCTTTTTTCTGATGTCCTGGAGTACCTTTC
 TCCTGGCGGAAGCCTGTGGGTTTACAGGTGTTGTAGCTGCTCTTTCTGTGGGATTACACAAGCTCATT
 CACATATAACAACCTGTGAGTAGAATCTCGAAGTGAAGCAAGCAGCTTTTTGAGGTGTTACACTTCCTG
 GCAGAGAACTTCTCTTTTCTACATGGCCTGGCACTGTTACCTTCCAGAAACACGTTTTTCAGCCCCA
 TTTTCATCATCGGTGCTTTTGTGCCATCTTCTGGGAAGAGCTGCACATATCTACCCACTCTCCTTCTT
 CCTCAACTTGGGCAGAAGGCATAAGATTGGCTGGAATTTTCAACACATGATGATGTTTTCAGGCCTGAGA
 GGAGCGATGGCATTGCACTGGCTATCCGTGACACGGCATCGATGCCCGCCAGATGATGTTCAACAACA
 CACTGCTCATTGTCTTCTCACTGTGTGGGTGATTGGTGGAGGCACGACCCCATGTTGTCATGGCTCAA
 TATAAGATACAGGTTCTGGAACATGATTGACAGA

ACGCGTACGCGGCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG208728 representing BC058750
 Red=Cloning site Green=Tags(s)


MEPSDAARPGPGRFRGLSPRLLLPLLPVLLGRGLRAGAAASSGEAAEDSSAMEELATEKEAEESHQRD
 SVSLLTFILLTLTILTIWLFKRRVRFLHETGLAMIYGLIVGVILRYGTPATSGHDKSLCTQEDRAFS
 TLLVNVSGKFFEYTLKGEISPGKINNVEQNDMLRKVTFDPEVFFNILLPPIIFHAGYSLKKRHFRLNLS
 ILAYAFGLTAVSCFIIGNLMYGVVKLMKIVGQLSDKFYTDCLFFGAIISATDPVTVLAIFNELHADVDL
 YALLFGESVLNDAVAIVLSSSIVAYQAPGLNTHAFDAAFFKSVGIFLGIFSGSFTMGAVTGVVTALVTK
 FTKLHCFPLLETALFFLMSWSTFLAEACGFTGVVAVLFCGITQAHYTYNNLSVESRSRSKQLFEVLHFL
 AENFIFSVMGLALFTFQKHVFSPIFIIGAFVAIFLGRAAHYPLSFFLNLRHKGWNFQHMMMFSGLR
 GAMAFALAIRDTASYARQMMFTTLLIVFFTVWVIGGGTTPMLSWLNIRYRFWNMIDR

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI


Cloning Scheme:



ACCN:	BC058750
ORF Size:	1644 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC058750.1
RefSeq Size:	2206 bp
RefSeq ORF:	1646 bp
Locus ID:	236727

Cytogenetics:	X A1.3
Gene Summary:	Mediates electroneutral exchange of protons for Na(+) and K(+) across endomembranes. May contribute to Golgi volume and cation homeostasis (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG208728